TECHNICAL INSTRUCTIONS & USER GUIDE

C100 CYLINDRICAL SHUTTLE LOCK

college park

TECHNOLOGY for the HUMAN RACE

FIGURE 1


FIGURE 2


FIGURE 3


FIGURE 4


FIGURE 5


FIGURE 6



FIGURE 7

PACKAGE CONTENTS (*Figure 1*)

- A. Lock Assembly
- B. Adjustable Length Push Pin
- C. (2) Ratchet Plunger Pins (*selected pin size*)
- D. (4) M6 x 25mm Fasteners (*not shown*)
- E. (1) Tube of CA Glue (*not shown*)

TOOLS RECOMMENDED

Loctite® 242

FABRICATION KIT (NOT INCLUDED) (*Figure 2*)

- A. Compression Nut
- B. Lamination Disk
- C. Fabrication Dummy
- D. SHSS 5/16-18 x 2.25" set screw
- E. Button Shield Dummy
- F. Adhesive Foam Pad
- G. Post Covers (*not shown*)

PRODUCT DESCRIPTION

The cylindrical shuttle lock includes a lock assembly with adjustable length push pin, two ratchet plunger pins, and four mounting fasteners.

INTENDED USE

The C100 Shuttle Lock is intended to secure a prosthetic liner to a socket as well as serve as an attachment point for endoskeletal components for lower limb prosthesis users.

INDICATIONS: CONTRAINDICATIONS:

- Lower limb amputations
- None known

TECHNICAL SPECIFICATIONS

Build Height	24 mm (0.95 in)
Assembly Weight	69 g (no pin)
Patient Weight Limit	136 kg (300 lbs)
Warranty	2 years

EN

PLUNGER PIN OPTIONS

Thread Type: Metric M10

Apply Loctite 242 and install into locking liner. Use a wrench to tighten.

PIN SIZES		
XS	5/8 in	16 mm
S	1 in	25 mm
STD	1.5 in	38 mm
L	2 in	51 mm
XL	2.5 in	64 mm

INSTRUCTIONS

PREPARE MODEL

1. Place the fabrication dummy on the distal end of the model, along the mid line and trace the circumference. (Figure 3A)
2. Create a flat surface on top of cast, to accommodate the size of the dummy.
3. Drill 1/4" (6 mm) hole in center of cast and embed set screw. Hex end of screw should protrude at least 1" (2.5 cm) outside of cast. (Figure 3B)
4. Continue on with fabrication method of choice: Thermoforming or Lamination

FABRICATION METHOD - THERMOFORMING

1. Add stockinet over model. Use string to tie off around set screw and remove excess material. (Figure 4A)
2. Thread fabrication dummy onto set screw, with correct button position. (Figure 4A)
3. Install button shield dummy and place adhesive foam pad over hex head hole. (Figure 4B)
4. Cover 4-hole posts with post covers. (Figure 4B)
5. Heat plastic using standard method. Then, thermoform over fabrication dummy. (Figure 4C)
Optional: A 2nd layer of plastic may be applied to the distal end for reinforcement.
Optional: The lamination disk may be used to create an even distal surface.
6. Sand distal end of the socket until the heads of the post covers are barely exposed and a flat surface has been created.
Caution: Inner and outer surfaces of distal socket must be flat and parallel in order to prevent lock from binding or breaking!
7. Expose push pin shield dummy and remove.
8. Extract model from socket.
9. Carefully sand around push pin hole area to create a smooth opening for push pin. (Figure 4D)
10. Continue to Lock Assembly.

FABRICATION METHOD - LAMINATION (FIGURE 5)

1. Apply resin deterrent to set screw threads. Then screw fabrication dummy onto it, with correct button position. (Figure 5A)
2. Apply inner PVA bag over dummy and model. Tie off PVA bag. (Figure 5B)
Note: Temporarily fitting the lamination disk and retaining nut will protect the PVA bag from the posts and set screw.
3. Create a seal around the dummy using PVA tape. (Figure 5B)
4. Trim excess PVA bag below push pin hole. Apply more tape to reinforce the seam. (Figure 5C)
5. Remove lamination disk and retaining nut. (Figure 5C)
6. Apply resin deterrent to the threads and hex head hole of button shield dummy and install. (Figure 5C)
7. Place adhesive foam pad on the button shield dummy. (Figure 5C)
8. Proceed with lay up materials. Ensure that posts and button shield dummy are exposed. (Figure 5D)
9. Place the lamination disc onto set screw, then tighten with compression nut, compressing disc against 4-hole posts. (Figure 5E)
10. Apply outer PVA bag and laminate socket. Ensure that the resin saturates all lay up material. (Figure 5F)
11. Once cured, remove excess resin and PVA bag.
12. Remove compression nut and lamination disk.
13. Expose push pin shield dummy and remove.
14. Remove model from socket.
15. Carefully sand the distal end of the socket to remove excess material. (Figure 5G)
16. Sand around push pin hole to create a smooth opening for push pin.

LOCK ASSEMBLY (FIGURE 6)

1. Insert lock assembly into bottom of socket, aligning push pin hole.
2. Align 4-hole pattern with socket and attach 4-hole adapter of choice.
3. Apply Loctite and torque all screws according to manufacturer's specifications.

PUSH PIN ASSEMBLY (FIGURE 7)

- A. Mark the desired length of the push pin stem and cut off excess material.
- B. Use CA glue to attach push pin cap to stem. Apply Loctite 242 and thread push pin into lock assembly.

 **WARNING**

- Apply Loctite® 242® on all mounting fastener threads.
- Use recommended torque settings to avoid malfunction.
- Do not expose this product to corrosive materials.
- Do not disassemble housing assembly or lock.
- Do not use lubricant in the locking mechanism, to prevent attraction of dirt and debris.

Failure to follow these guidelines may result in one or all of the following: noise, movement, damage, and catastrophic failure. For further questions, please contact College Park Technical Service.

USER GUIDE

Prosthetists: Instruct your patient on these steps and then give this page to the patient.

LOCK MAINTENANCE AND CARE

Inspect lock daily for performance and safety

Keep lock clean and dry to ensure smooth operation

Avoid moisture and humidity to prevent corrosion

Report any lock malfunction immediately to your prosthetist

WARRANTY INSPECTION / MAINTENANCE INFORMATION

College Park recommends that you schedule your patients for check-ups according to the Warranty Inspection schedule below.

High patient weight and/or impact level may require more frequent inspections. The C100 Cylindrical Shuttle Lock should be inspected after the first 30 days of use. We recommend you visually inspect the entire prosthesis for excessive wear and fatigue at each consultation.

WARRANTY INSPECTION SCHEDULE FOR C100 CYLINDRICAL SHUTTLE LOCK: 30 DAYS, THEN ANNUALLY.

TECHNICAL ASSISTANCE / EMERGENCY SERVICE 24-7-365

College Park's regular office hours are Monday through Friday, 8:30 am – 5:30 pm (EST). After hours, an emergency Technical Service number is available to contact a College Park representative.

LIABILITY

The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

CAUTION

College Park products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards are achieved only when College Park products are used with other recommended College Park components. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients.

CAUTION

If any problems occur with the use of this product, immediately contact your medical professional. The prosthetist and/or patient should report any serious incident* that has occurred in relation to the device to College Park Industries, Inc. and the competent authority of the Member State in which the prosthetist and/or patient is established.

*“Serious incident” is defined as any incident that directly or indirectly led, may have led, or might lead to any of the following: (a) the death of a patient, user, or other person, (b) the temporary or permanent serious deterioration of a patient’s, user’s, or other person’s state of health, (c) a serious public health threat.

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©2022 College Park Industries, Inc. All rights reserved.
College Park is a registered trademark.