AXIS[®]
college park

locking shoulder joint

technical instructions

التعليمات الفنية • Teknisk vejledning • Technische Anleitung • Τεχνικές οδηγίες • Instrucciones técnicas
Tekniset ohjeet • Instructions techniques • הוראות טכניות • Istruzioni tecniche • Technische instructies
Tekniske instruksjoner • Instrukcje techniczne • Instruções técnicas • Instruções técnicas • Технические
инструкции • Technické pokyny • Tekniska anvisningar • Teknik Talimatlar • 技术说明

FIGURE 1


FIGURE 3


FIGURE 2


FIGURE 4


FIGURE 5



FIGURE 6

LANGUAGES

EN-English	7
AR-Arabic	11
BG-Bulgarian	15
DA-Dansk	19
DE-Deutsch	23
EL-Elliniká	27
ES-Español	31

FI-Suomen Kieli	35
FR-Français	39
HR-Hrvatski	43
HU-Magyar.....	47
IT-Italiano	51
NL-Nederlands	55
NO-Norsk.....	59

PL-Polski.....	63
PT-Português.....	67
RO-Română.....	75
RU-Russkiy.....	79
SL-Slovenščina.....	83
SV-Svenska.....	87
TR-Türkçe.....	91
ZH-Zhōngwén.....	95

PACKAGE CONTENTS

(1) AXIS Shoulder Joint	(10) M3 0.5 x 20 mm PH FL screws	(10) M3 0.5 x 30mm PH FL screws
(1) 0.050" hex key	(1) SJ90 Mounting Ring	(1) 0.062" hex key
(1) 0.156" hex key	(1) 3mm hex key	(1) T10 Torx Key
(1) 0.125" Drill Bit	(1) SS Roll Pin (0.125" x 0.250")	

TOOLS REQUIRED

3mm drill bit	14mm drill bit
#1 Phillips head screwdriver	

ENDO ONLY

Torque wrench

This diagram (Figure 1) is to help familiarize you with the unique parts of the AXIS. These parts are referenced in the instructions and used when speaking with a technical service representative.

Key Components (Figure 1)

A. Exoskeletal Humeral Plate	B. Abduction Friction Adjustment	C. Yoke Plate	D. Mechanical Locking Mechanism	E. Electronic Locking Mechanism (Electric Model)
F. Size Mounting Ring	G. T10 Torx Humeral Plate Fasteners	H. Tube Clamp (Endo Model)	I. Endo Tube (Endo Model)	

AXIS Part Number	Mounting Type	Lock Actuator
SJ90 BK	Exoskeletal	Manual
SJ100 BK	Endoskeletal	Manual
SJ92 BK	Exoskeletal	Electronic
SJ102 BK	Endoskeletal	Electronic

PRODUCT DESCRIPTION

The Axis Locking Shoulder Joint is constructed with an exoskeletal or endoskeletal mounting option, and a mechanical or electronic lock option. It includes a mounting ring and ten mounting fasteners.

INTENDED USE

The Axis Shoulder Joint is a device intended to replace one or more functions of the biological human shoulder joint for upper limb amputees.

**INDICATIONS:**

- Upper limb amputations at the shoulder

**CONTRAINDICATIONS:**

- None known

TECHNICAL SPECIFICATIONS

Adduction	180°
Flexion-Extension	240° (25 locking positions)
Weight Capacity	23 kg (50 lbs) with prosthesis
Product Weight	290 g
Warranty	2 years

ENDOSKELETAL TUBE INSTRUCTIONS

1. Insert tube into tube clamp and secure by hand tightening the set screw. (Fig.2A)
2. Measure and cut to desired length.
3. Using the provided 0.125" drill bit, drill through the hole in the center of the tube clamp and through the tube itself. (Fig. 2B)
4. In the hole you've just created, fully insert the included SS roll pin. This will prevent the tube from pulling out. (Fig. 2C)
5. Torque the set screw of the tube clamp to 5 N-m.

MOUNTING INSTRUCTIONS (ALL MODELS)

1. Create room for mounting plate on plaster cast and attach
2. Pull check socket over cast and mounting plate
3. Clean off mounting plate for use in attaching shoulder to socket
4. Place mounting plate into recessed area of check socket
5. Drill holes for mounting screws
6. Attach shoulder joint to mounting plate using provided fasteners and torque to 1 N-m (9 in-lb)
7. Adjust abduction friction to desired level by turning the hex head screw at the top of the yoke clamp.
8. Humeral Plate can be removed to access wire routing channel. When re-installing plate, torque fasteners 1 N-m (9 in-lb).

MOUNTING INSTRUCTIONS (MECHANICAL)

Position the shoulder joint vertically so that the locking mechanism is pointed distally. If the locking mechanism is rotated too far clockwise or counterclockwise it will impede the Bowden cable of the remote lock or nudge kit. (Figure 3A)

Align the shoulder joint vertically when determining the mounting position. (Figure 3B)

Note: Avoid abduction angles that will cause contact with the yoke joint.

AXIS ELECTRONIC LOCK ACTUATOR

Note: Axis Electronic Shoulder Joints are compatible with the Steeper Espire Elbow Pro and Hybrid versions only.

Electronic Lock Actuator Overview

The Axis Shoulder Joint is available with an electronic lock actuator in both exoskeletal and endoskeletal versions. The shoulder is powered by the Espire Elbow battery and controlled via the input of choice. It requires an output cable to actuate the lock, which connects to the input board.

ROUTING CABLES THROUGH THE EXOSKELETAL SHOULDER JOINT (FIGURE 4)

In the exoskeletal version, both input and output wires are routed through a channel between the yoke plate and the humeral plate.

- (Figure 4A) To electrodes, switches, etc.
 - (Figure 4B) Lock output wire
 - (Figure 4C) To input board
1. A hole entering the socket must be provided for the Axis output wire.
 2. Input wires travel from the input control (electrodes, switches, etc) to the input connector board. The Axis output wire travels from the shoulder joint to the input connector board.
 3. Attach cables to input board.

Info: For more information on electrode fabrication, see **Espire Elbow Fabrication Instructions**. The Espire Elbow is manufactured by Steeper. Contact Steeper for more info:

ROUTING CABLES THROUGH THE ENDOSKELETAL SHOULDER JOINT (FIGURE 5)

- (Figure 5A) To electrodes, switches, etc.
- (Figure 5B) To input board
- (Figure 5B) Lock output wire

In the endoskeletal version, both input and output wires are routed through the center of the shoulder joint.

1. A hole entering the socket must be provided for the Axis output wire.
2. Run all cables from socket, through center of shoulder joint, down outside of pylon, to elbow input board.

3. Fasten cables to pylon with tape.
4. Insert cables through one of two openings on the endo adapter.
5. Attach cables to input board.

Info: For more information on endo adapter assembly, see **Espire Elbow Endo Adapter Kit Instruction Sheet**

Info: Use of a cosmetic cover is recommended to encase and protect cable wires.

CHANGING THE DEFAULT ROTATION STARTING POINT

The swing rotation can be changed for either a left-sided or right-sided amputation.

1. Remove fastener (Figure 6A)
2. Unlock AXIS, and rotate to desired position. (Figure 6B)
3. Re-Insert fastener into available slot. (Figure 6C)

TECH TIPS

- Use 30mm screws for attaching to a thick check socket and the 20mm for a thinner check socket. Recommended amount of screw length projecting out of the socket should be between .450-.550"
- Mounting surface must be flat and parallel in order for lock to operate properly.
- Do not remove the backing plate on the shoulder joint.
- The 1/16" hex-head screw, accessible through the center of the joint can be removed and relocated to change the position of the center stop. There are 6 possible locations, 60° apart.

- With the bent plate, the user can internally rotate 15° toward the body mid-line when the upper arm is straight ahead. If you prefer a straight line, then cut off the bent end of the Humeral Plate and mount the straight end to the Yoke.

WARNING

Load capacity may not exceed 50lbs (22.7 Kg) including prosthesis.

WARRANTY INSPECTION / MAINTENANCE INFORMATION

College Park recommends that you schedule your patients for check-ups according to the Warranty Inspection schedule below.

Inspection Guidelines

- Visual/Audible signs of wear, cracks, damage, etc.
- Ensure that the joint moves freely without obstruction.
- Ensure that the lock engages and disengages properly.

WARRANTY INSPECTION SCHEDULE FOR AXIS: ANNUALLY.

TECHNICAL ASSISTANCE / EMERGENCY SERVICE 24-7-365

College Park's regular office hours are Monday through Friday, 8:30 am – 5:30 pm (EST). After hours, an emergency Technical Service number is available to contact a College Park representative.

LIABILITY

The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer

CAUTION

College Park products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards are achieved only when College Park products are used with other recommended College Park components. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients.

CAUTION

If any problems occur with the use of this product, immediately contact your medical professional. The prosthetist and/or patient should report any serious incident* that has occurred in relation to the device to College Park Industries, Inc. and the competent authority of the Member State in which the prosthetist and/or patient is established.

*“Serious incident” is defined as any incident that directly or indirectly led, may have led, or might lead to any of the following; (a) the death of a patient, user, or other person, (b) the temporary or permanent serious deterioration of a patient's, user's, or other person's state of health, (c) a serious public health threat.

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©2023 College Park Industries, Inc. All rights reserved.