# COLLEGE PARK CAPITAL

## hydraulic knee

ركبة هيدروليكية • hydraulic knæ • hydraulisches Knie • υδραυλικό γόνατο  
rodilla hidráulica • hydraulinen polvi • Prothèses de genoux hydrauliques  
ברך הידראולית • ginocchio idraulico • hydraulische knie • hydraulisk kne  
hydrauliczne kolano • joelho hidráulico • joelho hidráulico  
гидравлический коленный сустав • hydraulický kolenný klb  
hydrauliskt knö • hidrolisk diz • 液压膝盖







# college park


TECHNOLOGY for the HUMAN RACE

## technical instructions


التعليمات الفنية • Teknisk vejledning • Technische Anleitung • Τεχνικές οδηγίες • Instrucciones técnicas  
Tekniset ohjeet • Instructions techniques • הוראות טכניות • Istruzioni tecniche • Technische instructies  
Tekniske instruksjoner • Instrukcje techniczne • Instruções técnicas • Instruções técnicas  
Τεχνικές ενσχυρισεις • Technické pokyny • Tekniska anvisningar • Teknik Talimatlar • 技术说明


**FIGURE 7**


**FIGURE 8**

**FIGURE 9**

**FIGURE 10**

# LANGUAGES

EN-English .....	7
AR-Arabic .....	14
BG-Bulgarian .....	21
DA-Dansk .....	28
DE-Deutsch .....	35
EL-Elliniká .....	42
ES-Español .....	49
FI-Suomen Kieli .....	56
FR-Français .....	63
HE-תִּירְבֵּעַ .....	70
HR-Hrvatski .....	78
HU-Magyar .....	84
IT-Italiano .....	91
NL-Nederlands .....	98
NO-Norsk .....	105
PL-Polski .....	112

PT-Português .....	119
PT-Português Brasil .....	126
RO-Română .....	133
RU-Russkiy .....	140
SK-Slovak .....	147
SV-Svenska .....	154
TR-Türkçe .....	161
UKR - українська .....	168
ZH-Zhōngwén .....	175

**PACKAGE CONTENTS**

(1) Capital Knee

**TOOLS RECOMMENDED**

(1) 4 mm Hex Key

(1) 2.5 mm Hex Key

*This diagram is to help familiarize you with the unique parts of the Capital Knee. These parts are referenced in the instructions and used when speaking with a technical service representative.*

**KEY COMPONENTS (FIGURE 1)****Both Models**

- A.** Pyramid (Proximal Attachment)
- B.** Upper Assembly
- C.** Knee Frame
- D.** Mode Adjustment Switch
- E.** Hydraulic Adjustment Cover
- F.** Knee Pad

**300 Model**

- G.** 34 mm Pylon Receiver (Distal Attachment)
- H.** 4 mm Clamp Screw (Torque 10 N-m)

**350 Model**

- I.** Pyramid (Distal Attachment)

**PRODUCT DESCRIPTION**

The Capital Knee is constructed with an integrated pyramid (proximal) and 34 mm pylon receiver or pyramid (distal).

**INTENDED USE**

The Capital Knee, intended for transfemoral amputees, is a prosthetic device designed to restore some function of an anatomical knee joint.

**INDICATIONS:**

Lower limb amputations above the knee

**CONTRAINDICATIONS:**

None known

## TECHNICAL SPECIFICATIONS

<b>Material</b>	Aluminum	<b>Assembly Weight</b>	980 g (34 mm Pylon Receiver) 990 g (Pyramid)
<b>Joint Type</b>	Single-Axis	<b>Patient Weight Limit</b>	150 kg (330 lbs)
<b>Flexion</b>	130°	<b>Warranty</b>	3 Years
<b>Build Height</b>	1.7 cm (0.67 in)	<b>Torque (Clamp Screw)</b>	10 N·m (7.4 ft-lbs)

### BUILD HEIGHT (FIGURE 2)

Distance from the axis of rotation to the proximal dome: 17mm (0.67 in)

### BENCH ALIGNMENT (TKA) (FIGURE 3)

- A. Determine the heel height
- B. Determine the socket flexion
- C. The alignment reference line bisects through the knee center

### KNEE CENTER REFERENCE (FIGURE 4)

Cross-hairs (⊕) indicate axis of rotation. Weight line can be anterior or posterior depending on desired overall alignment required.

### STATIC ALIGNMENT

With the alignment reference line through the knee center, plantarflex or dorsiflex the foot until the load line is balanced between 1/3 heel and 2/3 toe lever.

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Note: If the load line is too far anterior to knee center, it may become too difficult to initiate knee flexion. If the load line is too far posterior to knee center, it may cause premature knee flexion.

Caution: Using a 4-hole adapter for a long residual limb can limit knee flexion, because it may hit the knee frame before the full range of knee flexion is achieved. Consideration of appropriate components should be given during fabrication. After fabrication and before fitting, ensure there is not interference of lamination componentry and the knee frame.

FOR LONG RESIDUAL LIMBS, COLLEGE PARK RECOMMENDS • CPI 3-Prong Socket Adapter, 3PA RR T or CPI 4-Prong Socket Adapter, 4PA RR T

## ADJUSTMENTS

### 3-MODE FUNCTION (FIGURE 5)

The Capital Knee has three functional modes that can be adjusted by the end user. To change the mode, press down on the mode switch lock, then rotate the switch to the left, center, or right position. When lined up, the lock will click in place.

#### Lock Mode (Figure 6A)

When the switch is moved to the left, the knee provides high resistance to flexion. It will continue to extend until full extension is reached.

#### Normal Walking Mode (Figure 6B)

At the center position, the knee functions according to the hydraulic settings of swing flexion, swing extension, and stance flexion resistance.

#### Free Swing Mode (Figure 6C)

When the switch is moved to the right, stance flexion is bypassed. This allows free movement in both flexion and extension.

Caution: To prevent sudden collapse, users should use caution when applying weight to the knee during free swing mode.

### ACCESSING HYDRAULIC ADJUSTMENTS (FIGURE 7)

Remove the hydraulic adjustment cover before adjusting and replace when finished. (Figure 7A)

Note: Hydraulic valve adjustments can be made using a 4mm Allen wrench.

Note: The full range of valve adjustment is 180°. The effect of a dynamic adjustment can be felt with as little as 1/8 - 1/4 turn.

### DYNAMIC ADJUSTMENTS

For dynamic adjustments, have your patient begin ambulating at a comfortable cadence to fine-tune the flexion and extension resistance.

After the initial settings are established, it is important for the patient to demonstrate walking at variable speeds and on uneven surfaces (ramps, stairs, etc.) in order to fully adjust the knee. Perform a stand-to-sit action and/or stair/ramp descent to optimize stance flexion resistance.

**STANCE PHASE****STANCE FLEXION RESISTANCE (FIGURE 9C)**

**Factory setting = screw is at maximum resistance (most resistance to stance yielding) (Figure 8C)**

Controls flexion of the knee during stance phase of gait. This adjustment is most effective for controlling resistance from stand-to-sit, stair or ramp descent, and “stumble recovery.”

SYMPTOM	DESIRED RESULT	SCREW ADJUSTMENT	
Not Enough Flexion Resistance	Increase Stance Flexion Resistance	Turn <b>S</b> Clockwise	
Too Much Flexion Resistance	Decrease Stance Flexion Resistance	Turn <b>S</b> Counterclockwise	

**STANCE TRIGGER RESISTANCE (FIGURE 9D)**

**Factory setting = screw is at minimum resistance (easiest to trigger stance flexion) (Figure 10)**

This sets the amount of toe load required to trigger release into flexion. If the knee goes into flexion too easily at higher walking speeds, this adjustment may need to be increased.

SYMPTOM	DESIRED RESULT	SCREW ADJUSTMENT	
Swing Flexion Triggered Prematurely	Increase Trigger Resistance	Turn <b>T</b> Clockwise	
Swing Flexion Difficult To Trigger	Decrease Trigger Resistance	Turn <b>T</b> Counterclockwise	

**SWING PHASE**

**Factory setting = both screws at minimum resistance**

Begin by adjusting the flexion resistance, or heel rise. Then adjust the extension resistance to control terminal impact.

Caution: Flexion and extension must be possible at all settings.

**FLEXION RESISTANCE (FIGURE 9A)**

**Factory setting = screw is at minimum resistance (Figure 8A)**

Controls hydraulic resistance to knee flexion upon off-loading of the prosthesis. This adjustment can be used to limit heel rise.

SYMPTOM	DESIRED RESULT	SCREW ADJUSTMENT	
Flexion Too Fast or Excessive Heel Rise	Increase Resistance	Turn <b>F</b> Clockwise	
Flexion Too Slow or Heel Rise Not Sufficient	Decrease Resistance	Turn <b>F</b> Counterclockwise	

**EXTENSION RESISTANCE (FIGURE 9B)**

**Factory setting = screw is at minimum resistance (Figure 8B)**

Controls how the knee moves through extension during swing. This adjustment is most effective for reducing terminal impact.

SYMPTOM	DESIRED RESULT	SCREW ADJUSTMENT	
Extension Too Fast or Impact Too Abrupt	Increase Resistance	Turn <b>E</b> Clockwise	
Extension Too Slow or Not Sufficient	Decrease Resistance	Turn <b>E</b> Counterclockwise	

## ENVIRONMENTAL CONDITIONS

Operating Temperature: -7 to 55° C (20 - 130° F)

## USE IN WATER

The Capital Knee has been approved for use in fresh water.

- The knee should be switched into lock mode when using it in or near water.
- Use caution when walking on wet surfaces.
- After the knee encounters moisture, wipe it dry using a lint-free cloth.

## WARNING

- Flexion and extension must be possible at all settings.
- Large sockets must not contact the knee cylinder during flexion.
- Avoid pinching hazards! Do not place fingers near the flexing area of the knee.
- Patient concerns about the function should be reported to the prosthetist immediately, including but not limited to noise, sudden loss of function, etc.
- Do not disassemble the knee. Contact College Park to arrange a repair or replacement.
- Do not expose this product to corrosive materials, saltwater, or pH extremes.
- Contaminants such as dirt and the use of lubricants or powder may affect the function of the knee and lead to failure.
- Do not use compressed air to clean the knee as it can push dirt inside of the knee.

Failure to follow these technical instructions or use of this product outside the scope of its Limited Warranty may result in injury to the patient or damage to the product.

## WARRANTY INSPECTION AND MAINTENANCE INFORMATION

College Park recommends that you schedule your patients for check-ups per the warranty inspection schedule below.

High patient weight or activity level may require more frequent inspections. We recommend you visually inspect the following applicable parts for excessive wear and fatigue at each warranty inspection.

- Knee Assembly, Hydraulic Cylinder, Knee Pad

Warranty inspection schedule for the Capital: Six months, then annually.

**ISO 10328 - P7 - 150 kg\*** 

*\* Body mass limit not to be exceeded!*

*For specific conditions and limitations of use, see intended use section of manufacturer's written instructions.*

## COMPLIANCE

This device has been tested according to ISO 10328 standard to three million load cycles. Depending on patient activity, this may correspond to 3-5 years of use.

## CAPITAL KNEE

The load level specification is ISO10328 – P7 – 150 kg

## TECHNICAL ASSISTANCE / EMERGENCY SERVICE 24-7-365

College Park's regular office hours are Monday through Friday, 8:30 am – 5:30 pm (EST). After hours, an emergency Technical Service number is available to contact a College Park representative.

## LIABILITY

The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer

## CAUTION

If the product is exposed to conditions outside of the recommended environmental conditions inspect for signs of changes or loss of functionality during use which may include poor response, lack of swing flexion initiation, abnormal noises, etc. If there are noticeable changes or loss of function discontinue use of the product.

## CAUTION

College Park products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards are achieved only when College Park products are used with other recommended College Park components. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients.

## CAUTION

If any problems occur with the use of this product, immediately contact your medical professional. The prosthetist and/or patient should report any serious incident\* that has occurred in relation to the device to College Park Industries, Inc. and the competent authority of the Member State in which the prosthetist and/or patient is established.

\*'Serious incident' is defined as any incident that directly or indirectly led, may have led, or might lead to any of the following; (a) the death of a patient, user, or other person, (b) the temporary or permanent serious deterioration of a patient's, user's, or other person's state of health, (c) a serious public health threat.

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©2025 College Park Industries, Inc. All rights reserved.