ODYSSEY iQ

microprocessor foot

college park

TECHNOLOGY *for the* HUMAN RACE

DOMESTIC
technical instructions

technische anweisungen • instrucciones técnicas
instructions techniques • istruzioni tecniche

FIGURE 1


FIGURE 2


FIGURE 3


FIGURE 4


FIGURE 5


FIGURE 6



FIGURE 7

LANGUAGES

EN-English	6
DE-Deutsch	29
ES-Español	52
FR-Français	75
IT-Italiano	98

This document is intended for use by certified and licensed prosthetic clinicians only.

For user approved operations, see the Odyssey iQ Microprocessor Foot User Guide.

Federal Law restricts that the Odyssey iQ microprocessor foot is to be purchased, configured, and fit only by a board-certified prosthetist, licensed by the state in which they practice. This device is intended for use in accordance with the information contained in this document. Instruct the patient on proper use of this device before transferring device to patient.

Thank you for purchasing the Odyssey iQ microprocessor Foot from College Park Industries. The following pages will cover everything you need to know about this system from setup to operation. If you have any questions, concerns, or comments, please contact our Technical Service team at 800 728 7950 (US/Canada), (+1) 586 294 7950 (International).

In the following document you will find information on everything from initial setup to maintenance and care of the Odyssey iQ microprocessor foot. Read these instructions carefully and educate the end user on all functions of this product before final delivery.

The Odyssey iQ microprocessor foot is fully assembled and undergoes electronic testing verification before it is shipped.

What's in the box?

Hardware

Odyssey iQ Microprocessor Foot
CPI Sock
Foot Shell
Battery Charger w/ Adapter

Instruction Manuals

Odyssey iQ Microprocessor Foot Technical Instructions
Odyssey iQ Microprocessor Foot Quick Setup Guide
Odyssey iQ Microprocessor Foot User Manual
Stride Studio App for Odyssey iQ Instructions

Tools Recommended:

(1) 4mm Hex Key
(1) Foot Horn

These diagrams are to help familiarize you with the unique parts of the Odyssey iQ foot. These parts are referenced in the instructions and used when speaking with a technical service representative.

KEY COMPONENTS (FIGURE 1)

A. Integrated Pyramid

D. Heel Spring

G. Battery

J. Serial Number

B. Ankle-Housing Assembly

E. Permanent Wedge

H. Power Button / LED Indicator

• CPI Sock (not shown)

C. Toe Spring

F. Rounded Heel Wedge (Optional)

I. Charging Receptacle / Dust Cover

• Foot Shell (not shown)

DRAIN HOLE

A pre-cut hole is located on the bottom of the foot shell to allow a draining path for water (*Figure 2*).

FUNCTION

The Odyssey iQ foot is a hydraulically dampened microprocessor-controlled foot. A system of integrated sensors is used to determine when to make changes to plantarflexion and dorsiflexion resistances. Functions include:

FUNCTIONS

- Level ground walking
- Changes in walking speed
- Ramp ascent/descent
- Low/critical battery mode

FEATURES

- Waterproof in fresh water (IP67)

USER SELECTED MODES

- Stance iQ mode
- Walking mode
- Ankle lock
- Custom shoe profiles
- Bilateral foot linking

PRODUCT DESCRIPTION

This prosthetic foot device is constructed with an integrated pyramid, ankle-housing assembly, two composite springs, heel wedge, and battery. The toe spring is secured to the housing and heel spring with fasteners.

INTENDED USE

The Odyssey iQ is a prosthetic foot designed to replace one or more functions of the biologic human foot.

INDICATIONS:

- Lower limb amputations

CONTRAINDICATIONS:

- Patients exceeding the 125 kg / 275 lbs weight limit.
- Unusual activities, including extreme sports.
- Use in corrosive environments, including salt or chlorinated water.

PROTECTIVE COVER ON DOME

Remove the protective cover on dome after alignment is completed and before patient leaves clinic.

TECHNICAL SPECIFICATIONS

STRUCTURAL SPECIFICATIONS

Foot Size	Weight Limit	Build Height	Foot Weight*
23-25 cm	275 lbs / 125 kg	3.5 in / 8.8 cm	1,027 g
26-30 cm		3.6 in / 9.0 cm	

*26cm foot w/shell

GAIT MATCHING GUIDELINES

The gait match determines the firmness of the foot based on the user's specifications (foot size, patient weight, and impact level).

FIRMNESS CATEGORIES

Refer to the chart below to determine the correct firmness category.

Note: Incorrect category selection may result in poor device function. Contact College Park Technical Service if you have questions about category selection.

FIRMNESS CATEGORY CHART

WEIGHT LBS	0-140	141-180	181-220	221-275
WEIGHT KG	0-63	64-81	82-100	101-125
SIZE CM	23-30			
LOW IMPACT	1	2	3	4
MODERATE IMPACT	2	3	4	5

DEVICE OPERATION

Battery	Lithium polymer, 3.7V, 980 mAh
Time to full charge	4 hours (time to 50% charge = 2 hours)
Charging Cable	5V, 3A DC (Use only provided power supply to charge the battery)

WIRELESS

Connection	Bluetooth 5.0/Bluetooth Low Energy (BLE)	Maximum Speed	24 Mbps
Operating Frequency	2.402 - 2.480 GHz	Effective Radiated Power	12.00 dBm

ENVIRONMENTAL USE CONDITIONS

Charging	32F to 104F (0C to +40C)
Operating	14F to 122F (-10C to +50C), Pressure: 70-106 kPa, Humidity: 0%-100% relative humidity
Storage & Transport	-4F to 140F (-20C to +60C) *, Pressure: 70-106 kPa, Humidity: 0%-90% relative humidity

Note: If storing device above or below operating temperature, allow the device to return to within operating temperature range before use. The device must be brought up to the operating temperature, allow the device to sit for 15 minutes.

IP RATING

IP67	Dust tight and protected against freshwater at a maximum depth of 1 meter for 30 minutes.
-------------	---

SAFETY

Only start up the product in accordance with the information contained in the supplied documents.

Failure to follow these technical instructions or use of this product outside the scope of its Limited Warranty may result in injury to the patient or damage to the product.

LEGEND OF SYMBOLS

	Note: Possible technical damage.	 CAUTION	Caution: Possible risk of accident or injury.
	Info: Basic information regarding this product.	 WARNING	Warning: Possible risk of severe accident or injury.

SAFETY INSTRUCTIONS

	Info: Use on airplanes Airlines may not permit the use of this device on their aircraft. Check with the airline before travelling to ensure this device is allowed for use on the plane.
	Info: Disposal These products may not be disposed of with household waste in some jurisdictions. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.
 CAUTION	Caution: Risk of pinching where the ankle joint bends Ensure that fingers and other body parts are not in this area when bending the ankle joint.
 CAUTION	Caution: Unsupervised use It is not recommended for children to operate this device without the supervision of an adult. Use extreme caution around small children and household pets.
 CAUTION	Caution: Risk of accident while operating a vehicle A lower extremity amputee's ability to drive a vehicle is determined on a case-by-case basis. Factors include the type of fitting (amputation level, unilateral or bilateral, residual limb conditions, design of the prosthesis) and the amputee's abilities. All persons are required to observe their country's national and state driving laws when operating vehicles. For insurance purposes, drivers should have their driving ability examined and approved by an authorized test center.

Caution: Improper use

Any type of excessive strain, overload or improper use may lead to faulty control or malfunction of the Odyssey iQ foot, resulting in a risk of injury. The Odyssey iQ foot was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports with excessive strain and/or shocks, or extreme sports. Careful handling of the prosthesis and its components not only increases their service life but, above all, ensures your personal safety! Should the prosthesis be subjected to unusual stresses (such as a fall), immediately contact a certified prosthetist and have the prosthesis inspected for any damage.

Caution: Consequences of product deterioration

Wear and tear on system components can lead to malfunction of the Odyssey iQ foot, resulting in a risk of injury. Follow the specified service intervals. The service life of this device is 3 years. Patient concerns about the function should be reported to the prosthetist immediately, including but not limited to noise, sudden loss of function, etc.

Caution: Mechanical overloading

External mechanical influences or loads, such as impacts and vibration, can lead to faulty control or malfunction of the Odyssey iQ foot and result in a risk of injury. The Odyssey iQ foot should not be subjected to mechanical vibrations or impacts.

Caution: Manipulation of system components

Independent changes and/or modifications to system components may lead to faulty control or malfunction of the Odyssey iQ foot, resulting in a risk of injury. No modifications of Odyssey iQ foot except those described in this information document are authorized. The Odyssey iQ foot and damaged components may only be opened or repaired by certified College Park technicians.

Caution: Battery damage

Damage can occur to the battery when dropping, knocking, crushing, vibrating, or puncturing. Avoid damaging lithium batteries and devices. Always inspect for signs of damage, such as hissing, leaking, cracking/ bulging, and smoking before use. If any of these signs are present, immediately place away from flammable materials. Remove the device from service and contact your prosthetist. Do not charge.

Caution: Penetration of dirt and moisture into charge port

The penetration of dirt or moisture into the charge port may lead to faulty control or malfunction of the Odyssey iQ foot and result in a risk of injury. Dirt/debris should be removed from the Odyssey iQ charge port. Ensure it is dry before charging.

Caution: Water and humidity

The foot can be operated in freshwater but cannot be used in extreme situations such as diving or jumping into water. It is protected against freshwater at a maximum depth of 1 meter for 30 minutes. The foot is not resistant to corrosion and should not come into contact with corrosive materials, saltwater, chlorinated water, or pH extremes..

CAUTION

- Use caution when walking on wet surfaces.
- Under certain circumstances, use in wet environments may increase the risk of slip and fall.
- The foot should be switched into ankle lock mode when using it in water or near wet surfaces.
- After contact with water, ensure that the water has drained out of the system and that the drainage port has not become blocked (Figure 2: drain hole in shell)
- After the foot encounters moisture, wipe it dry using a lint-free cloth.
- Make sure that the charge port is clean and dry before charging.
- If the foot encounters salt or chlorinated water, it should be rinsed immediately and dried off
- If the device malfunctions after exposure to water, contact College Park.

Caution: Thermal overloading

Extreme temperature conditions can lead to faulty control or malfunction of the Odyssey iQ foot and result in a risk of injury. Avoid areas outside the specified operating temperature range. The operating temperature range must be between **-10 °C and 50 °C (14 °F and 122 °F)**.

CAUTION

Caution: Critically low battery status or loss of power

The foot will go into safety mode when in a critically low state. Safety mode means that the foot goes to the set plantarflexion and dorsiflexion resistances and will no longer compensate for terrain or speed changes and won't be able to enter Stance iQ mode.

CAUTION

Caution: Magnetic interference

The Odyssey iQ foot and connected components can malfunction when near high-tension power lines, transmitters, transformers, or other sources of strong electromagnetic radiation (such as security systems for goods in department stores). This can result in a risk of injury.

CAUTION

Caution: Too close to HF communication devices (e.g. mobile phones, Bluetooth devices, WIFI devices)

The device may be susceptible to electromagnetic interference from portable and mobile RF communications devices such as mobile (cellular) telephones or other equipment, even if that other equipment complies with CISPR EMISSION requirements. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Odyssey iQ foot, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result. The Odyssey iQ foot is suitable for use in any environment except where immersion in fluids other than freshwater is possible, or where exposure to highly electrical and/or magnetic fields can occur (e.g. electrical transformers, high power radio/TV transmitters, RF surgical equipment, CT and MRI scanners).

CAUTION

**Caution: EM disturbances**

Do not use the Odyssey iQ foot near active HF SURGICAL EQUIPMENT and the RF shielded room of an EM SYSTEM for magnetic resonance imaging, where the intensity of EM DISTURBANCES is high. High levels of Electro Magnetic disturbances may cause system to stop functioning properly, either not responding to input signal(s) or no movement of the joints.

**Caution: Operating the product near active implanted systems**

When operating the product, there is a risk of temporary influences of active implantable systems (e.g. pacemakers, defibrillators etc.) because of electromagnetic interference of the product.

- **When operating the product in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.**
- **Make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.**

**Warning: Using with other equipment**

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

**Warning: Use only specified equipment**

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

**Warning: Operation of equipment in hospitals**

The Odyssey iQ foot was designed for use in residential environments (home, restaurants, etc.) not hospitals or industrial areas. If device is used in environments such as hospitals or industrial areas, user might have to relocate to operate device appropriate so that is away from other HF radio devices.

**Warning: Use in an oxygen rich environment**

The device is not designed for use in presence of flammable gases or oxygen rich environments.

THE BATTERY

The Odyssey iQ foot contains an internal lithium polymer battery. This battery supplies 3.7 V 980 mAh. It is advised to charge the battery daily. For most users, the battery will last up to 16 hours of continuous use, depending on the condition of the battery and the frequency of use.

The battery is shipped with a partial charge (up to 30%). We recommend charging the battery to 100% before use.

LED INDICATOR

The Odyssey iQ foot is equipped with a state-of-charge indicator. This indicator tells the user how much life is left in their current battery. The Odyssey iQ foot must be powered on to use this feature.

Press the power button for 1-second to activate the LED indicator light (*Figure 1-H*). The number of blinks indicates the state of the charge. Charge level can also be checked via the Stride Studio app.

FEATURE

POWER BUTTON

Description

Press and hold for 1 second to indicate the battery status.

FOOT LED INDICATIONS – BATTERY

COLOR

Solid Green

4 flashes

3 flashes

3 flashes

1 flash

INDICATOR

BATTERY STATUS

Fully charged, 91-100%

75 - 90%

50 - 74%

25 - 49%

1 - 24%

CRITICALLY LOW BATTERY STATUS

The foot will go into safety mode when in a critically low state. Safety mode means that the foot goes to the set plantarflexion and dorsiflexion resistances and will no longer compensate for terrain or speed changes and won't be able to enter Stance iQ mode.

BATTERY CHARGER

Odyssey iQ feet are supplied with a charger for the lithium polymer battery. The charger is recommended for daily use and will ensure that the battery will receive a full charge and provide maximum running time. There are four power adapter options (US, UK, EU, or AU) to match the needs of different regions. Use only the provided power supply to charge the lithium polymer battery.

CHARGING THE BATTERY

CHARGING THE ODYSSEY iQ FOOT BATTERY

1. Only charge the foot on a flat, level surface and in a cool, dry location away from sources of heat and moisture. Plug the wall adapter into the main AC supply (Figure 3-B).
2. To connect to the front of the ankle, open the dust cover and insert the connector.

Note: The connector is directionally specific (Figure 3-A).

3. The foot will audibly indicate when it begins charging or is unplugged.

Recharge time from empty is approximately 4 hours.

POWER BUTTON

The power button is located on the top of the ankle, near the anterior opening of the foot shell (Figure 1-H). To turn the device on or off, press and hold the power button for 3 seconds. When the device is powered on or off, a green light on the LED indicator will display for about 3 seconds.

FEATURE

POWER BUTTON

Description

Turns foot on/off. Initiates Bluetooth pairing when prompted in app.

POWER BUTTON

ON/OFF Press and hold the power button for 3 seconds.

PAIRING ODYSSEY iQ FOOT TO MOBILE DEVICE

BLUETOOTH PAIRING

The foot must be paired to the Stride Studio app when it is used for the first time. Launch the app, select "Connect New Device", and a list of available devices will be displayed. Select the serial number of the foot from the menu. A Bluetooth pairing request will pop up. Press the "Pair" button. The app will display an "Unauthorized" pop up. Press the power button on the front of the ankle for 1 second to complete the pairing process.

RECOMMENDED ALIGNMENT

ENDOSKELETAL MOUNTING

Use only high quality proximal endoskeletal components.

ASSEMBLY AND DISASSEMBLY (FOR SOCK REPLACEMENT)

Use a foot horn to don and doff the foot shell. Remove the CPI sock and replace as needed. Any further disassembly or modification of components will void the warranty.

Note: To avoid damage to the foot or battery, use care when donning or doffing the foot shell.

HYDRAULIC RANGE

The Odyssey iQ foot has 12° of hydraulic motion. The foot is designed to provide 2.5° hydraulic dorsiflexion from the neutral standing position (*Figure 4*).

The foot was developed utilizing a dynamic carbon fiber base. As a result, a typical user will experience an additional range of dynamic motion during ambulation.

Note: Excessive angular adjustment will affect the hydraulic range of the foot. After making an alignment change, ensure that the user retains 2.5° of hydraulic dorsiflexion.

SETUP WIZARD

If this is the first time setting up a new Odyssey iQ foot, the Setup Wizard will automatically launch. The Setup Wizard is designed to guide the prosthetist through the required steps of the setup process.

Order of Operation:

- Linking
- Bench & Static Alignment
- Calibrate Shoe Height
- Dynamic Adjustments
- Calibrate Walking Speed

LINKING

The Linking process is set up by the prosthetist, so that prosthetists and bilateral users can easily switch between the left foot and right foot during app use.

Select whether the setup is for a unilateral or bilateral foot. If the foot is unilateral, the app will skip the linking process and advance to the next section. If the foot is bilateral, follow the app instructions.

BENCH ALIGNMENT

The Odyssey iQ was designed with spring-assisted dorsiflexion, which dorsiflexes the ankle 4 degrees when unweighted.

1. The heel rise is 3/8" (10 mm). The load line divides the foot at 1/3 heel lever and 2/3 toe lever (*Figure 5A*).
2. Install the mounting adapter so that it aligns with the reference line on the ankle dome (*Figure 5B*).
3. Once bench alignment is completed, don the shoe. Pulling the ankle back 4 degrees will align it in a neutral position (*Figure 6*). The ankle will return to a dorsiflexed position if no weight is applied.

STATIC ALIGNMENT

Begin with the hydraulic valves at the factory settings (plantarflexion = 50%, dorsiflexion = 45%). Have the user stand comfortably and evaluate the heel-toe balance of the foot (*Figure 7*).

They will sense the hydraulic movement of the ankle but should not feel like they are falling forward or backward. Use alignment to position the foot at the point where they feel the most balanced.

SYMPTOM	ALIGNMENT CHANGE
Falling Backward	Shift foot posterior relative to the socket
Falling Forward	Shift foot anterior relative to the socket

CALIBRATE SHOE HEIGHT

This section calibrates the angle for shoes with heel heights between 0-0.80 in (0-2.0 cm).

To calibrate the shoe height, users should begin in a standing position, wearing their preferred shoes. Follow the app instructions to complete the process.

Note: Attempting to calibrate shoes outside of the specified range may result in an error message.

DYNAMIC ALIGNMENT

Have the user begin by walking on level ground, to evaluate the heel-toe resistance and gait timing. Using the Stride Studio App, adjust for plantarflexion resistance first, then dorsiflexion.

Finalize dynamic alignment by observing the user walking on an inclining-declining surface (ramp). Make further adjustments to the resistance valves as necessary.

Plantarflexion resistance affects the user's gait from Heel Strike to Foot Flat.

Dorsiflexion resistance affects the user's gait through Midstance, as the body travels over the foot.

DYNAMIC ADJUSTMENTS

DESIRED RESULT	ADJUSTMENT	COMPONENT CHANGE
Firmer Heel Response	Increase plantarflexion resistance	Install Rounded Heel Wedge
Softer Heel Response	Decrease plantarflexion resistance	Remove Rounded Heel Wedge
Firmer Toe Response	Increase dorsiflexion resistance	NONE
Softer Toe Response	Decrease dorsiflexion resistance	NONE

CALIBRATE WALKING SPEED

There is a calibration process to capture the user's normal walking speed. This is to improve accuracy in foot performance. Foot calibration should be performed after alignment is complete and all resistance adjustments have been made. To calibrate walking speed, follow the app instructions.

Note: Excessive amounts of external rotation may limit sensitivity to certain functions of the foot.

ADDITIONAL CONSIDERATIONS

Have the user practice standing up from a seated position in order to acclimate to the motion of the ankle. Use caution when driving. Make sure the user is comfortable with the motion of the ankle if using the Odyssey iQ as their driving foot or turn on ankle lock mode.

APP OVERVIEW

DOWNLOAD STRIDE STUDIO APP

College Park's Stride Studio App is available for both iOS & Android devices. Download Stride Studio from Apple App Store or Google Play Store.

MINIMUM OPERATING SYSTEMS:

- Apple iOS 15
- Android Version 8 (Oreo)

Select Stride Studio icon on the device to launch the app.

DASHBOARD SCREEN

SYMBOL	NAME	DESCRIPTION	SYMBOL	NAME	DESCRIPTION
	Setup	Alignment recommendations.		Battery Status	Indicates battery charge level of Odyssey iQ.
	Adjustments	Adjustments to Odyssey iQ foot resistance and mode selection.		Device Warnings	List of error codes.
	Shoes	Allows prosthetist and user to set custom shoe profiles.		Alerts Settings	Turn non-critical alert sounds on/off.
	Activity	Activity log: download and print reports.	OD/Q000001	Serial Number	Displays the serial number of the foot.
	Information	Technical instructions and online learning resources.		Link Foot	Allows bilateral users to switch easily between left and right feet.

USER SELECTED MODES

STANCE iQ MODE

Stance iQ mode is a function where the foot detects that the user is standing still and enters high resistance.

If Stance iQ mode is enabled, the foot will lock plantar and dorsiflexion for added support, after 8 seconds of standing still. Once the user begins walking again, it will revert to normal function allowing for plantar and dorsiflexion.

WALKING MODE

In this mode, Stance iQ is disabled and the foot will not enter high resistance.

ANKLE LOCK MODE

Enters the foot into high resistance until deactivated.

CUSTOM SHOE PROFILES

Allows prosthetist and user to program additional shoe profiles. Calibrates the angle for shoes with heel heights between 0-0.80 in (0-2.0 cm).

BILATERAL FOOT LINKING

Set up by the prosthetist, so that prosthetists and bilateral users can easily switch connection to the left or right foot during app use.

Info: For information on use with mobile devices, refer to the **Stride Studio App Instructions**.

TROUBLESHOOTING

Caution: The Odyssey iQ foot should never be serviced while connected to the end user. Ensure that the device is disconnected and powered off before any service or maintenance is performed. This device should never be serviced while in use. Never let children use this handle this device unsupervised. Take caution when using this device around pets that may cause damage to the device.

The Odyssey iQ foot features an LED indicator located on the top of the ankle, near the anterior opening of the foot shell. This LED is used to display battery related information.

The chart below outlines what each of the different LED patterns means. When the device is turned on the LED will blink briefly. Once the device is on, battery life can be checked by pressing the power button for 1 second.

ODYSSEY iQ FOOT ALL LED INDICATIONS

COLOR	INDICATOR	BATTERY STATUS
Solid Green		Fully charged, 91-100%
4 flashes		75 - 90%
3 flashes		50 - 74%
3 flashes		25 - 49%
1 flash		1 - 24%

If the Odyssey iQ is not functioning as expected, try the following:

- Turn the system off, wait several seconds and power back on.
- Make sure the battery has sufficient charge. If the battery is too low, charge the battery.

WARRANTY INSPECTION AND MAINTENANCE INFORMATION

CLEANING INSTRUCTIONS

Clean the Odyssey iQ foot by wiping it with a soft, lint-free cloth. Never apply harsh chemicals.

USE IN WATER

Clean the Odyssey iQ foot by wiping it with a soft, lint-free cloth. Never apply harsh chemicals.

- The Odyssey iQ foot has been approved for use in freshwater at a maximum depth of 1 meter for 30 minutes.
- Use caution when walking on wet surfaces.
- The foot should be switched into ankle lock mode when using it in water or near wet surfaces.
- Always dry foot thoroughly after exposure to moisture. Rinse and dry if exposed to salt water or chlorinated water.

COSMETIC FINISHING

College Park does not recommend cosmetically finishing the prosthesis if it will be exposed to moisture.

WARRANTY INSPECTION AND MAINTENANCE INFORMATION

College Park recommends that you schedule your patients for check-ups per the warranty inspection schedule below.

High patient weight or activity level may require more frequent inspections. We recommend you visually inspect the following applicable parts for excessive wear and fatigue at each warranty inspection.

- Hydraulic Assembly
- Composites and Adapters
- Battery
- CPI Sock
- Foot Shell
- Wedge

WARRANTY INSPECTION SCHEDULE FOR THE ODYSSEY iQ: SIX MONTHS, THEN ANNUALLY.

TECHNICAL ASSISTANCE / EMERGENCY SERVICE 24-7-365

College Park's regular office hours are Monday through Friday, 8:30 am – 5:30 pm (EST). After hours, an emergency Technical Service number is available to contact a College Park representative

LIABILITY

The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

CAUTION

College Park products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards are achieved only when College Park products are used with other recommended College Park components. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients.

CAUTION

If any problems occur with the use of this product, immediately contact your medical professional. The prosthetist and/or patient should report any serious incident* that has occurred in relation to the device to College Park Industries, Inc., and the competent authority of the Member State in which the prosthetist and/or patient is established.

*‘Serious incident’ is defined as any incident that directly or indirectly led, may have led, or might lead to any of the following: (a) the death of a patient, user, or other person, (b) the temporary or permanent serious deterioration of a patient's, user's, or other person's state of health, (c) a serious public health threat.

COMPLIANCE

This device has been tested according to ISO 10328 standard to two million load cycles. The expected life duration under normal usage is 2 million steps. Depending on the amputee's activity this corresponds to a duration of use of two to three years.

ISO 10328 - LABEL

WEIGHT LIMIT (KG)	LABEL TEXT
125	ISO 10328-P6-125 kg

ISO 10328 - “P” - “m”kg*) ⚠

***) Body mass limit not to be exceeded!**

For specific conditions and limitations of use, see intended use section of manufacturer's written instructions.

SYMBOLS

FEATURE	DESCRIPTION	FEATURE	DESCRIPTION
	Direct current	IC: 23763- RSL10SIP	IC ID of BT Module
	Caution		FCC Mark 21 CFR Part 15
	See operating instructions/consult instructions for use		Non-ionizing radiation
	Type BF Applied Part		CE Mark
	Refer to instruction manual/booklet		Device contains electronic components and/or batteries that should not be disposed of in regular waste
IP67	Dust tight and protected against freshwater at a maximum depth of 1 meter for 30 minutes.		On/Off
	Keep away from water/keep dry		Recyclable Lithium Polymer Battery
	Temperature limitation	Lithium	
	Medical device manufacturer		Medical Device
	Lot number		Complies with Australia Radio communications requirements
	Serial number		Note: Possible technical damage
	Humidity Rating		Info: Basic information regarding this product
	Pressure Rating		Caution: Possible risk of accident or injury
FCC ID: 2APD9- RSL10SIP	FCC ID of BT Module		Warning: Possible risk of severe accident or injury
		CAUTION	
		WARNING	

COMPLIANCE

STANDARD	EDITION	DESCRIPTION
ISO 10328	Structural testing of lower-limb prostheses-requirements and test methods	2016
ANSI/AAMI 60601-1	Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance	AAMI ES60601-1: 2005 & A1:2012 & A2:2021
IEC 60601-1-2	Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests	2014 Ed.4+A1
IEC 60601-1-6	Medical Electrical Equipment - Part 1-6: General Requirements for Basic Safety and Essential Performance - Collateral Standard: Usability	2010 Ed. 3+A1;A2
IEC 62366-1	Medical Devices - Part 1: Application of Usability Engineering to Medical Devices	IEC 62366-1 ED. 1.1 B:2020
IEC 60601-1-11	Medical Electrical Equipment - Part 1-11: General Requirements for Basic Safety and Essential Performance - Collateral Standard: Requirements for Medical Electrical Equipment and Medical Electrical Systems Used in the Home Healthcare Environment	2015 Ed.2+A1
IEC 62304	Medical Device Software - Software Life Cycle Processes	IEC 62304 ED. 1.1 B:2015
FCC Part 15	Radio Frequency	
IEC 62133	Secondary Cells And Batteries Containing Alkaline Or Other Non-Acid Electrolytes - Safety Requirements For Portable Sealed Secondary Lithium Cells, And For Batteries Made From Them, For Use In Portable Applications - Part 2: Lithium Systems	IEC 62133-2 ED. 1.1 B:2021
IEC 60529	Degrees of Protection Provided by Enclosures (IP Code)	2013 Ed.2.2

EMC COMPLIANCE – SPECIFIC MITIGATIONS

The Odyssey iQ foot was tested to the listed standards at the appropriate levels for Home Health Care Equipment below to ensure safety of the product regarding immunity and emissions. All devices maintained their performance during and after the tests were completed.

This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference. (2) This device must accept any interference received, including interference that may cause undesired operation.

This device complies with Industry Canada license-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

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PHENOMENON & STANDARD	TEST LEVEL	REMARKS
Radiated Emissions CISPR11 ed5.0 (with A1:2010), CISPR 11 ed6.1 (2015 +A1:2016)	Group 1, Class B	<i>The Odyssey iQ foot uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.</i>
Electro-Static Discharge Immunity Test IEC 61000-4-2 ed2.0 (2008-12)	Contact ± 8 kV Air ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV	<i>The Odyssey iQ foot is attached to the patients socket which is design by a certified prosthetist.</i>
Radiated, Radio-Frequency, Electromagnetic Immunity IEC 61000-4-3 ed3.0 (with A1:2007+A2:2010)	10 V/m 80 MHz – 2,7 GHz 80% AM at 1 kHz	<i>Portable and mobile RF communications equipment should be used no closer to any part of the Odyssey iQ foot, including cables, than the recommended in the technical manual. Device does not need to be operated in a shielded environment.</i>
Power Frequency Magnetic Field Immunity Test IEC 61000-4-8 ed2.0 (2009-09)	30 A/m, 50 Hz or 60 Hz	<i>The Icon Knee should not be operated closer than 15 cm of sources of The Odyssey iQ foot should not be operated closer than 15 cm of sources of power frequency magnetic field.</i>

TEST SPECIFICATIONS FOR ENCLOSURE PORT IMMUNITY TO RF WIRELESS COMMUNICATIONS EQUIPMENT

TEST FREQUENCY (MHZ)	BAND (MHZ)	SERVICE	MODULATION	MAXIMUM POWER (W)	DISTANCE (M)	IMMUNITY TEST LEVEL (V/M)
385	380-390	TETRA 400	Pulse modulation 18 Hz	1,8	0,3	27
450	430-470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0,3	28
710	704-787	LTE Band 13, 17	Pulse modulation 217 Hz	0,2	0,3	9
745						
780						
810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0,3	28
870						
930						
1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1. 3. 4, 25; UMTS	Pulse modulation 217 Hz	2	0,3	28
1845						
1970						
2450	2400-2570	Bluetooth WLAN, 802.11 b/g/n, RFIO 2450, LTE Band 7	Pulse modulation 217 Hz	2	0,3	28
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0,2	0,3	9
5500						
5785						

TEST SPECIFICATIONS FOR ENCLOSURE PORT IMMUNITY TO PROXIMITY MAGNETIC FIELDS

TEST FREQUENCY	MODULATION	IMMUNITY TEST LEVEL (A/M)
30 kHz	CW	8
134,2 kHz	Pulse modulation 2,1 kHz	65
13,56 MHz	Pulse modulation 50 kHz	7,5

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©2024 College Park Industries, Inc. All rights reserved.