college park
TECHNOLOGY for the HUMAN RACE

technical instructions

التعليمات الفنية • Teknisk vejledning • Technische Anleitung
Τεχνικές οδηγίες • Instrucciones técnicas • Tekniset ohjeet
Instructions techniques • הוראות טכניות • Istruzioni tecniche
Technische instructies • Tekniske instruksjoner • Instrukcje techniczne
Instruções técnicas • Instruções técnicas • Технические инструкции •
Technické pokyny • Tekniska anvisningar • Teknik Talimatlar • 技术说明

FIGURE 1


FIGURE 2


FIGURE 3



FIGURE 4

LANGUAGES

EN-English	4
AR-Arabic	7
BG-Bulgarian	10
DA-Dansk	13
DE-Deutsch	16
ES-Español	19
FI-Suomen Kieli	22
FR-Français	25
HR-Hrvatski	28
HU-Magyar	31
IT-Italiano	34
JA-Nihongo	37
KO-Kugo	40
NL-Nederlands	43
NO-Norsk	46
PL-Polski	49
PT-Português	52
PT-BR-Português Brasil	55
RO-Română	58
RU-Russkiy	61
SL-Slovenščina	64
SV-Svenska	67
TR-Türkçe	70
UK-Ukrayins'ka	73
ZH-Zhōngwén	76

Package Contents

- (1) Odyssey K3 Foot (1) Foot Shell
- (1) CPI Sock (1) Odyssey K3 Wedge Kit
- (1) Cosmetic Attachment Plate (CAP)

Tools Recommended

- (1) 3mm Hex Key
- (1) 4mm Hex Key
- (1) Foot Horn

This diagram is to help familiarize you with the unique parts of the Odyssey K3. These parts are referenced in the instructions and used when speaking with a technical service representative.

Key Components (Figure 1)

- A. Integrated Pyramid
 - B. Ankle-Housing Assembly
 - C. Toe Spring
 - D. Heel Spring
 - E. Permanent Wedge
 - F. Rounded Heel Wedge (Optional)
 - G. Nut Plate
- CPI Sock (not shown)
 - Foot Shell (not shown)

PRODUCT DESCRIPTION

This prosthetic foot device is constructed with an integrated pyramid, ankle-housing assembly, two composite springs, and heel wedge. The toe spring is secured to the housing and heel spring with fasteners.

INTENDED USE

The Odyssey K3 is a prosthetic foot designed to replace one or more functions of the biologic human foot.

INDICATIONS:

- Lower limb amputations

CONTRAINDICATIONS:

- None known

PROTECTIVE COVER ON DOME

Remove the protective cover on dome after alignment is completed and before patient leaves clinic.

TECHNICAL SPECIFICATIONS

FOOT SIZE	WEIGHT LIMIT	BUILD HEIGHT	FOOT WEIGHT*
21-25 cm	275 lbs / 125 kg	2.95 in / 7.5 cm	749 g
26-30 cm		3.04 in / 7.7 cm	

*26cm foot w/shell

GAIT MATCHING® GUIDELINES

The gait match determines the firmness of the foot based on the user's specifications (foot size, patient weight, and impact level).

FIRMNESS CATEGORIES

Refer to the chart below to determine the correct firmness category.

Note: Incorrect category selection may result in poor device function. Contact College Park Technical Service if you have questions about category selection.

WEIGHT LBS	0-140	141-180	181-220	221-275
WEIGHT KG	0-63	64-81	82-100	101-125
SIZE CM	21-30			
LOW IMPACT	1	2	3	4
MODERATE IMPACT	2	3	4	5

MOUNTING

Use only high quality endoskeletal components.

ASSEMBLY AND DISASSEMBLY (FOR SOCK REPLACEMENT)

Use the Foot Horn to don and doff the foot shell. Remove the CPI Sock and replace as needed. Any further disassembly or modification of components will void the warranty.

HYDRAULIC RANGE (Figure 2)

The Odyssey K3 foot has 12° of hydraulic motion. The foot is designed to provide 3° hydraulic dorsiflexion from the neutral standing position.

The foot was developed utilizing a dynamic carbon fiber base. As a result, a typical user will experience an additional range of dynamic motion during ambulation.

Note: Excessive angular adjustment will affect the hydraulic range of the foot. After making an alignment change, ensure that the user retains 3° of hydraulic dorsiflexion.

STATIC ALIGNMENT (Figure 3)

For optimal function, balance the patient's weight evenly between the heel and toe.

- The Odyssey K3 was designed with a 3/8" (10 mm) heel rise.
- The load line divides the foot at 1/3 heel lever and 2/3 toe lever.

Using a 3mm hex key for adjustment, the hydraulic valves should be set at minimum resistance. Have the user stand comfortably and evaluate the heel-toe balance of the foot.

They will sense the hydraulic movement of the ankle but should not feel like they are falling forward or backward. Use alignment to position the foot at the point where they feel the most balanced.

SYMPTOM	ALIGNMENT CHANGE
Falling Backward	Shift foot posterior relative to the socket
Falling Forward	Shift foot anterior relative to the socket

DYNAMIC ADJUSTMENTS (Figure 4)

Have the user begin by walking on level ground, to evaluate the heel-toe resistance and gait timing. Using a 3mm hex key, adjust for planterflexion resistance first, then dorsiflexion.

Finalize dynamic alignment by observing the user walking on an inclining-declining surface (ramp). Make further adjustments to the resistance valves as necessary.

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Planterflexion resistance affects the user's gait from Heel Strike to Foot Flat.

Dorsiflexion resistance affects the user's gait through Midstance, as the body travels over the foot.

DESIRED RESULT	VALVE ADJUSTMENT	COMPONENT CHANGE
Firmer Planterflexion	Turn P-valve clockwise (more resistance)	Install Rounded Heel Wedge
Softer Planterflexion	Turn P-valve counterclockwise (less resistance)	Remove Rounded Heel Wedge
Firmer Dorsiflexion	Turn D-valve clockwise (more resistance)	NONE
Softer Dorsiflexion	Turn D-valve counterclockwise (less resistance)	NONE

ADDITIONAL CONSIDERATIONS

Have the user practice standing up from a seated position in order to acclimate to the motion of the ankle. Use caution when driving. Make sure the user is comfortable with the motion of the ankle if using the Odyssey K3 as their driving foot.

WARNING

- Do not expose this product to corrosive materials, salt water or pH extremes.
- Contaminants such as dirt and the use of lubricants or powder may effect the function of the CPI Sock and lead to noise.
- Failure to follow these technical instructions or use of this product outside the scope of its Limited Warranty may result in injury to the patient or damage to the product.

RESIDUAL RISK STATEMENT

NOTICE OF RESIDUAL RISK

During fitting process, ensure that CPI sock does not become pinched between foot and endoskeletal componentry.

WARRANTY INSPECTION /MAINTENANCE INFORMATION

College Park recommends that you schedule your patients for check-ups according to the warranty inspection schedule below.

High patient weight and/or impact level may require more frequent inspections. We recommend you visually inspect the following applicable parts for excessive wear and fatigue at each warranty inspection.

- Hydraulic Assembly
- Composites and Adapters
- CPI Sock
- Foot Shell
- Wedge

WARRANTY INSPECTION SCHEDULE FOR ODYSSEY K3 : 6 MONTHS, THEN ANNUALLY.

TECHNICAL ASSISTANCE / EMERGENCY SERVICE 24-7-365

College Park's regular office hours are Monday through Friday, 8:30 am – 5:30 pm (EST). After hours, an emergency Technical Service number is available to contact a College Park representative.

LIABILITY

The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

CAUTION

College Park products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards are achieved only when College Park products are used with other recommended College Park components. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients.

CAUTION

If any problems occur with the use of this product, immediately contact your medical professional. The prosthetist and/or patient should report any serious incident* that has occurred in relation to the device to College Park Industries, Inc. and the competent authority of the Member State in which the prosthetist and/or patient is established.

*"Serious incident" is defined as any incident that directly or indirectly led, may have led, or might lead to any of the following; (a) the death of a patient, user, or other person, (b) the temporary or permanent serious deterioration of a patient's, user's, or other person's state of health, (c) a serious public health threat.

COMPLIANCE

This device has been tested according to ISO 10328 standard to two million load cycles. Depending on patient activity this may correspond to 2-3 years of use.

ISO 10328 - LABEL

FOOT SIZE	WEIGHT LIMIT (KG)	LABEL TEXT
21-30 cm	125	ISO 10328-P6-125 kg

ISO 10328 - "P" - "m"kg*)

***) Body mass limit not to be exceeded!**

For specific conditions and limitations of use, see intended use section of manufacturer's written instructions.

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©2023 College Park Industries, Inc. All rights reserved.