pneumatic knee

رُكْبَة هوائية • pneumatisk knæprotese • Pneumatisches Knie
πνευματική άρθρωση γόνατος • rodilla neumática
pneumaattinen polvi • genou pneumatique • ברך פנאומטית
ginocchio pneumatico • pneumatische knie • pneumatisk kne •
kolano pneumatyczne • joelho pneumático • Joelho pneumático
пневматический коленный модуль • pneumatický kolenný klb
pneumatisk knä • pnömatik diz • 气动膝盖

technical instructions

التعليمات الفنية • Teknisk vejledning • Technische Anleitung
Τεχνικές οδηγίες • Instrucciones técnicas • Tekniset ohjeet
Instructions techniques • הוראות טכניות • Istruzioni tecniche
技術説明書 • 기술 지침서 • Technische instructies
Tekniske instruksjoner • Instrukcje techniczne • Instruções técnicas
Instruções técnicas • Технические инструкции • Technické pokyny
Tekniska anvisningar • Teknik Talimatlar • 技术说明

FIGURE 1


FIGURE 2


FIGURE 3





FIGURE 4

PACKAGE CONTENTS

(1) Victor Knee

TOOLS RECOMMENDED

(1) 4 mm Hex Key

EN

This diagram (Figure 1) is to help familiarize you with the unique parts of the Victor Knee. These parts are referenced in the instructions and used when speaking with a technical service representative.

Key Components, KNEE (Figure 1)

A. Pyramid

B. Upper Assembly

C. Lower Assembly

D. 30mm Pylon Receiver

E. 4mm Clamp Screw (Torque 10 N-m)

PRODUCT DESCRIPTION

The Victor Knee is constructed with an integrated pyramid (proximal) and 30 mm pylon receiver (distal).

INTENDED USE

The Victor Knee, intended for transfemoral amputees, is a prosthetic device designed to restore some function of an anatomical knee joint.

INDICATIONS:

- Lower limb amputations

CONTRAINDICATIONS:

- None known

TECHNICAL SPECIFICATIONS

Material	Aluminum
Joint Type	Single-Axis
Flexion	145°
Build Height	2.9 cm (1.1 in)
Assembly Weight	710 g
Patient Weight Limit	125 kg (275 lbs)
Torque (Clamp Screw)	10 N-m (7.4 ft-lbs)

CLEARANCE (Figure 2)

A	Overall height	8.0 in (20.4 cm)
B	Dome to knee center	1.1 in (2.9 cm)
C	Dome to tube end contact	6.3 in (16.0 cm)

BENCH ALIGNMENT (Figure 3)

- Determine the heel height
- Determine the socket flexion
- The load line bisects through the pylon

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ADJUSTMENTS (Figure 4)

EN All adjustments can be made using a 4mm Allen wrench.

STANCE PHASE

Load Dependence (Figure 4A)

Controls the brake sensitivity. When load is applied, the knee will not bend until the load is displaced. Turning the load adjustment clockwise will increase the amount of load required to initiate braking.

Factory setting = Screw is at minimum (most brake sensitivity).

Note: The factory setting for load dependence should generally not require adjustment. However, it may be necessary in cases such as patients with higher body weights, more aggressive gait, or alignment differences.

Stance Adjustment (Figure 4B)

Controls the amount of stance flexion. Adjustments can be made to alter the amount of stance flexion during loading response.

Factory setting = Turn screw clockwise to the point where resistance is felt, then back out about 1/8 turn (counterclockwise).

Readjustment after continued use – The screw should generally not require adjustment; however, play may develop over time. Turn screw 1/8 clockwise to reduce play.

Note: Stance adjustment works together with load dependence to control the brake function. If decreasing the stance adjustment, it may be necessary to also increase the load dependence adjustment.

⚠ Caution: Do not over tighten. May lead to brake sticking which can cause unstable gait.

SWING PHASE

Flexion Resistance (Figure 4C)

Controls the rate of flexion.

Extension Resistance (Figure 4D)

Controls the rate of extension.

Factory setting = both screws at minimum (counterclockwise).

Begin by adjusting the flexion resistance or heel rise first, then the extension resistance to control terminal impact.

Note: After adjusting the swing phase resistance, it may be necessary to readjust the brake.


⚠ Caution: Flexion and extension must be possible at all settings.

DYNAMIC ADJUSTMENTS

Note: The effect of a dynamic adjustment can be felt with as little as 1/8 - 1/4 turn.

STANCE PHASE

Load Dependence

SYMPTOM	DESIRED RESULT	SCREW ADJUSTMENT	
Brake function too strong; Swing too difficult to initiate	Increase Load Dependence (brake less sensitive)	Turn L clockwise	
Brake function not sufficient; Knee too unstable during stance	Decrease Load Dependence (brake more sensitive)	Turn L counterclockwise	

Stance Adjustment

EN

SYMPTOM	DESIRED RESULT	SCREW ADJUSTMENT	
Not enough flexion	Increase Stance Flexion	Turn S counterclockwise	
Too much flexion	Decrease Stance Flexion	Turn S clockwise	

SWING PHASE

Flexion Resistance

SYMPTOM	DESIRED RESULT	SCREW ADJUSTMENT	
Flexion too fast or Excessive heel rise	Increase Resistance (close valve)	Turn F clockwise	
Flexion too slow or Heel rise not sufficient	Decrease Resistance (open valve)	Turn F counterclockwise	

Extension Resistance

SYMPTOM	DESIRED RESULT	SCREW ADJUSTMENT	
Extension too fast or Impact too abrupt	Increase Resistance (close valve)	Turn E clockwise	
Extension too slow or Not sufficient	Decrease Resistance (open valve)	Turn E counterclockwise	

WARNING

- Flexion and extension must be possible at all settings.
- Avoid pinching hazards! Do not place fingers near the flexing area of the knee.
- Patient concerns about the function should be reported to the prosthetist immediately, including but not limited to: noise, sudden loss of function, brake release sticking, etc.
- Do not disassemble the knee. Contact College Park to arrange a repair or replacement.
- Do not expose this product to moisture, i.e. fresh water, corrosive materials, salt water or pH extremes. If the knee encounters moisture, wipe it dry using a lint-free cloth.
- Contaminants such as dirt and the use of lubricants or powder may affect the function of the knee brake and lead to failure.
- Do not use compressed air to clean the knee as it can push dirt inside of the knee.
- Failure to follow these technical instructions or use of this product outside the scope of its Limited Warranty may result in injury to the patient or damage to the product.

WARRANTY INSPECTION AND MAINTENANCE INFORMATION

College Park recommends that you schedule your patients for check-ups per the warranty inspection schedule below.

High patient weight or activity level may require more frequent inspections. We recommend you visually inspect the following applicable parts for excessive wear and fatigue at each warranty inspection.

- Knee Assembly, Extension Stop Bumpers

Warranty inspection schedule for College Park Victor: Six months, then annually.

TECHNICAL ASSISTANCE / EMERGENCY SERVICE 24-7-365

College Park's regular office hours are Monday through Friday, 8:30 am – 5:30 pm (EST). After hours, an emergency Technical Service number is available to contact a College Park representative.

LIABILITY

The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer

CAUTION

College Park products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards are achieved only when College Park products are used with other recommended College Park components. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients.

CAUTION

If any problems occur with the use of this product, immediately contact your medical professional. The prosthetist and/or patient should report any serious incident* that has occurred in relation to the device to College Park Industries, Inc. and the competent authority of the Member State in which the prosthetist and/or patient is established.

*'Serious incident' is defined as any incident that directly or indirectly led, may have led, or might lead to any of the following: (a) the death of a patient, user, or other person, (b) the temporary or permanent serious deterioration of a patient's, user's, or other person's state of health, (c) a serious public health threat.

COMPLIANCE

This device has been tested according to ISO 10328 standard to three million load cycles. Depending on patient activity this may correspond to 3-5 years of use.

ISO 10328 - LABEL

WEIGHT LIMIT (KG)

125

LABEL TEXT

ISO 10328-P6-125 kg

ISO 10328 - "P" - "m"kg*)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use, see intended use section of manufacturer's written instructions.

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©2023 College Park Industries, Inc. All rights reserved.