



Partial Hand
SOLUTIONS, LLC



USER GUIDE

M-FINGER
PARTIAL M-FINGER
M-THUMB

⚠ CAUTIONS/WARNINGS

- Do not expose this product to corrosive or caustic materials, salt water, or pH and temperature extremes.
- Do not disassemble.



Partial Hand Solutions
86 Lady Slipper Lane
Southington, CT USA 06489
+1 (860) 538-5532



EMERGO EUROPE
Prinsessegracht 20,
2514 AP
The Hague, Netherlands



M-FINGERS



INSPECTION/MAINTENANCE

- Users should clean the inside of the prosthesis and silicone liner with fragrant-free, anti-bacterial wipes daily.
- Check attachment screws to ensure they are tight.
- Evaluate cabling to make sure it is not fraying.
- Inspect grip tip for wear and replace as needed
- Visual/Audible signs of wear, cracks, damage,etc.
- Ensure that the joint moves freely without obstruction.

POSITIONING THE FINGERS

- Flex wrist to bend fingers
- Extend wrist to extend fingers

ADDITIONAL CONSIDERATIONS/TECH TIPS

- New users should start with a wearing schedule and ease in to a full day of use.
- If using a silicone socket, apply a small amount of water or Cal Stat®lubricant to limb to dampen skin. Push limb into prosthesis.
- If a strap is used on the forearm section, secure in place.
- Adjust BOA reel to the ideal cable tension. Flex wrist to bend fingers, extend wrist to extend fingers.
- To remove the prosthesis, remove the strap, loosen the BOA reel tension and remove limb from prosthesis.

PARTIAL M-FINGERS



INSPECTION/MAINTENANCE

- Users should clean the inside of the prosthesis and silicone liner with fragrant-free, anti-bacterial wipes daily.
- Visual/Audible signs of wear, cracks, damage, etc.
- Ensure that the joint moves freely without obstruction.
- Inspect grip tip for wear and replace as needed.

POSITIONING THE PARTIAL FINGER

- Flex knuckle (MCP joint) to flex finger
- Extend knuckle (MCP joint) to extend finger

ADDITIONAL CONSIDERATIONS/TECH TIPS

- New users should start with a wearing schedule and ease in to a full day of use.
- Apply wristband to hand. Apply small amount of water or Cal Stat® lubricant to dampen finger/fingers and push each finger into the prosthetic sockets.
- To remove partial fingers from limb, pull on socket to remove. Remove wrist band.

M-THUMBS



INSPECTION/MAINTENANCE

- Users should clean the inside of the prosthesis and silicone liner with fragrant-free, anti-bacterial wipes daily.
- Visual/Audible signs of wear, cracks, damage,etc.
- Ensure that the joint moves freely without obstruction.
- Inspect grip tip for wear and replace as needed.

POSITIONING THE THUMB

- Push on the back of the M-Thumb to flex the digit by bending the joint.
- Push on the inside of the M-Thumb to straighten the joint.
- The M-Thumb can be rotated 360 degrees

ADDITIONAL CONSIDERATIONS/TECH TIPS

- New users should start with a wearing schedule and ease in to a full day of use.
- If using a silicone socket, apply a small amount of water or Cal Stat® lubricant to limb to dampen skin. Push limb into prosthesis.