

# AXIS<sup>®</sup> SHOULDER JOINT

Formerly the LTI Locking Shoulder Joint, the new AXIS<sup>®</sup> has been re-engineered by College Park. With a friction hinge, the patented shoulder joint is operated with simple switches that simplifies harnessing by eliminating the need for gross body movements. The AXIS is able to be integrated with all upper limb control systems with a special channel to conceal and protect power cables. Unlike any product on the market, the AXIS provides better task control, a natural range of motion, and ease of use for the patient's daily life.



Endo model

SWING	LOCKING POSITIONS	ABDUCTION	WEIGHT LIMIT
240°	25	180°	50 lbs
L-CODE*	WARRANTY	MOUNTING	WEIGHT
L6646	2 years	Endo Exo	290g

†The listing of codes with these products should not be construed as a guarantee for coverage or payment. Ultimate responsibility for the coding of services/products rests with the individual practitioner.



240° Range of motion



**PART NUMBER  
EXO MODEL  
SJ90 BK**

**PART NUMBER  
ENDO MODEL  
SJ100 BK**




#### INCLUDED

AXIS Shoulder joint,  
mounting ring and  
hardware.

## SWITCHES

	PART NUMBER	DESCRIPTION	ADDITIONAL INFORMATION
	SJ43	Lock Release Kit	
	SJ45	Nudge Kit	

## ADDITIONAL HARDWARE

	PART NUMBER	DESCRIPTION	ADDITIONAL INFORMATION
	SJ60	3 Spoke Mounting Plate, Black	
	SJ63	Straight Mounting Plate, Aluminum	
	SJ53	Boston Arm Endo Adapter Kit	



# ORDER FORM - AXIS SHOULDER JOINT

NORTH AMERICA

PAGE 1

**CUSTOMER #** **P.O.#** **DATE**

## PROSTHETIST INFORMATION

### BILLING

FACILITY/ATTN:

ADDRESS

CITY STATE/ZIP

PHONE FAX

CARRIER\* ☐ UPS ☐ OTHER

\*ORDERS SHIP UPS GROUND ON THE FOLLOWING DAY UNLESS SPECIFIED.

### SHIPPING (LEAVE BLANK IF SAME AS BILLING)

FACILITY/ATTN:

ADDRESS

CITY STATE/ZIP

PHONE FAX

DATE REQUIRED TIME

EMAIL

## PATIENT INFORMATION

PROSTHETIST NAME

REQUISITIONER

PATIENT ID / NAME

NOTES

**PART NUMBER**  
**EXO MODEL**  
**SJ90 BK**

**PART NUMBER**  
**ENDO MODEL**  
**SJ100 BK**

### INCLUDED

AXIS Shoulder joint,  
mounting plate and  
hardware.

## ADDITIONAL PARTS

PART NUMBER	DESCRIPTION
SJ60	3 Spoke Mounting Plate, Black
SJ63	Straight Mounting Plate, Aluminum
SJ43	Lock Release Kit
SJ45	Nudge Kit
SJ53	Boston Arm Endo Adapter Kit

PART NUMBER	QTY



27955 College Park Dr, Warren, MI 48088, USA

EMAIL COMPLETED ORDER FORM TO: [orders@college-park.com](mailto:orders@college-park.com)  
TOLL FREE 800.728.7950 FAX 586.294.0067  
ON THE WEB [www.college-park.com](http://www.college-park.com)

573 FRM OF SJ 170127