

C100 *cylindrical shuttle lock*



The College Park C100 cylindrical shuttle lock is made from a high-strength durable nylon and used as a traditional lock system for lower limb prosthesis users. Retrofittable with similar shuttle locks, the C100 is lightweight, secure, and the push pin can easily be cut to the appropriate length for the prosthesis. The C100 lock is rated to 300 lbs with a 2-year warranty and comes with two ratchet plunger pins.

C100 PART NUMBER

C100 SLK

PIN SIZE

SEE CHART

PIN SIZES	
XS	5/8 in
S	1 in
STD	1.5 in
L	2 in
XL	2.5 in

INCLUDED

C100 Shuttle Lock, 2 ratchet plunger pins (indicate size), 4 M6 x 25mm fasteners and push pin kit.

TECHNICAL SPECIFICATIONS

FUNCTION	Suspension	WARRANTY	2 years
PIN SIZES	XS, S, STD, L, XL	WEIGHT	88 g (with pin)
MOUNTING	4-hole adapter	BUILD HEIGHT	0.95 in
WEIGHT LIMIT	300 lbs	L-CODE*	L5672

*The listing of codes with these products should not be construed as a guarantee for coverage or payment. Ultimate responsibility for the coding of services/products rests with the individual practitioner.

ADDITIONAL HARDWARE	PART NUMBER	DESCRIPTION	ADDITIONAL INFORMATION
	RP XS	Shuttle Lock Ratchet Pin 2 pk, 5/8"	Additional pair of ratchet pins
	RP S	Shuttle Lock Ratchet Pin 2 pk, 1"	
	RP STD	Shuttle Lock Ratchet Pin 2 pk, 1.5"	
	RP L	Shuttle Lock Ratchet Pin 2 pk, 2"	

ADDITIONAL HARDWARE	PART NUMBER	DESCRIPTION	ADDITIONAL INFORMATION
	RP XL	Shuttle Lock Ratchet Pin 2 pk, 2.5"	Additional pair of ratchet pins
	C100 SLK FAB	C100 Shuttle Lock Fabrication Kit	Kit for assistance with building socket for C100 Shuttle Lock.
	SL PPK	Shuttle Lock Push Pin Kit	Replacement push pin for C100 Shuttle Lock
	SL FRK	Shuttle Lock Fabrication Refill Kit	
	4HFK 25	M6 x 25mm Fasteners, Qty 4	