

EASY MODE SELECTION | SIMPLE ADJUSTMENTS | WATERPROOF IN FRESH WATER



The single-axis Capital hydraulic knee is easy to use for both patients and clinicians. Featuring smooth and responsive hydraulics during both swing and stance phase, this knee was designed to effortlessly accommodate the variable gait of the user. Durable, sleek, and waterproof in fresh water, the Capital is handmade in America with care. This knee also features easy mode selection, stumble assist technology, and standard 4mm adjustments for simple set up. It is a strong, reliable solution for moderate to high activity users.

CAPITAL PART NUMBER

CHKA

COLOR*	DISTAL ATTACHMENT**
G, J, R OR W	300 OR 350



*(G)rey, (J)et Black, (R)ed, or (W)hite.



** (300) 34 mm Pylon Receiver, (350) Pyramid




TECHNICAL SPECIFICATIONS


MATERIAL	Aluminum
DISTAL ATTACHMENT	34mm Pylon Receiver or Pyramid
FLEXION RANGE	130°
WEIGHT	1010 g (Pylon Receiver) 1020 g (Pyramid)
WEIGHT LIMIT	330 lbs
WARRANTY	2 years
DOMES TO KNEE CENTER	0.67 in
OVERALL HEIGHT	9.0 in
L-CODES*	L5828, L5845, L5850, L5925
WATER CAPABILITY  Waterproof in Fresh Water	

*The listing of codes with these products should not be construed as a guarantee for coverage or payment. Ultimate responsibility for the coding of services/products rests with the individual practitioner.



RECOMMENDED ADAPTERS (34 MM PYLON RECEIVER)		
	PART NUMBER	DESCRIPTION
	TA S250 34	250 mm Low-Profile Tube Adapter, 34 mm diameter, Stainless Steel
	TA S450 34	450 mm Low-Profile Tube Adapter, 34 mm diameter, Stainless Steel

RECOMMENDED ADAPTERS (PYRAMID)		
	PART NUMBER	DESCRIPTION
	TA T250 30	250 mm Low-Profile Tube Adapter, 30 mm diameter, Titanium
	TA T430 30	430 mm Low-Profile Tube Adapter, 30 mm diameter, Titanium
	TCT 30	Tube Clamp, 30 mm diameter, Titanium

HARDWARE/TOOLS		
	PART NUMBER	DESCRIPTION
	CP AJ	Knee Adjustment Cover
	CP KTCS 20	Knee Tube Clamp Screw, M5 x 20mm
	CP HK 4	Hex Key, 4mm

REPLACEMENT PARTS		
	PART NUMBER	DESCRIPTION
	CP CRK	Cylinder Replacement Kit
	CP KPJ	Knee Pad Replacement Kit, Jet Black
	CP KPW	Knee Pad Replacement Kit, White
	CP KPG	Knee Pad Replacement Kit, Gray
	CP KPR	Knee Pad Replacement Kit, Red
	CP LCJ	Knee Lower Cover Replacement Kit, Jet Black
	CP LCW	Knee Lower Cover Replacement Kit, White
	CP LCG	Knee Lower Cover Replacement Kit, Gray
	CP LCR	Knee Lower Cover Replacement Kit, Red