

**INTELLIWEAVE COMPOSITE | ACTIVE HEEL | 300 LB WEIGHT LIMIT**


The Celsus foot brings College Park's innovative Integrated Spring Technology (iST) into a low-impact design. It features an active heel and excellent ground contact to provide smooth, stable transitions. Durable yet lightweight, the Celsus was created to provide anatomical function for low-impact patients.

**CELSUS PART NUMBER**

<b>CS</b>	SIDE	SHELL COLOR*	<b>EN</b>	SIZE	<b>D</b>	FIRMNESS CATEGORY
	L OR R	C, T, OR B		21-30		1-4

**INCLUDED**




Endo foot includes: sandal toe foot shell and CPI sock

\*(C)aucasian, (T)an or (B)rown. Example part number: CS LCEN 24 D3

**CELSUS FIRMNESS CATEGORIES**

PATIENT WEIGHT	0-140lbs	141-180	181-220	221-300
SIZE	21-30cm			25-30cm
LOW IMPACT	1	2	3	4

**TECHNICAL SPECIFICATIONS**

<b>FUNCTION</b>	K2 iST
<b>SIZES</b>	21-30 cm
<b>MOUNTING</b>	Endo (IP)
<b>WEIGHT LIMIT</b>	220 lbs (21-24 cm) 300 lbs (25-30 cm)
<b>WARRANTY</b>	2 years
<b>BUILD HEIGHT</b>	2.2 in (21-24 cm) 2.5 in (25-26 cm) 2.6 in (27-30 cm)
<b>WEIGHT*</b>	526 g
<b>L-CODE**</b>	<b>IMPACT LEVEL</b>
L5972	
<b>WATER CAPABILITY</b>	<b>ENVIROSHELL***</b>
 Weatherproof	

**ENVIROHELL PART NUMBER**


<b>CS S</b>	SIDE	SHELL COLOR**	SIZE
	L OR R	C, T, OR B	21-30




\*\* (C)aucasian, (T)an or (B)rown




sandal toe

\*26 cm with Enviroshell. \*\*The listing of codes with these products should not be construed as a guarantee for coverage or payment. Ultimate responsibility for the coding of services/products rests with the individual practitioner. \*\*\*Caucasian, Tan, or Brown.



HARDWARE	PART NUMBER	DESCRIPTION
	CP SSET	Super Set Screws (8pc) M8 x 14mm
	TFH 800	Foot Horn, Stainless Steel
	CP HK 4	Hex Key, 4mm

SHELLTREAD	PART NUMBER (LEFT/RIGHT)	DESCRIPTION	ADDITIONAL INFORMATION
	CS TLN x CS TRN x	Shelltread (left or right), Neutral (indicate size)	Increases durability and prevents premature wear of foot shell.

SOCKS	PART NUMBER	DESCRIPTION	ADDITIONAL INFORMATION
	CPI SS SM	CPI Sock, Small, 21-25cm	Highly durable sock that acts as a noise barrier and protects the foot from debris.
	CPI SS LG	CPI Sock, Large, 26-30cm	
	CPI SSK SM	CPI Sock, Small, 21-25cm, 6-pack	
	CPI SSK LG	CPI Sock, Large, 26-30cm, 6-pack	