

CARBON COMPOSITE | HIGH WEIGHT LIMIT | ULTRA LOW PROFILE



The Horizon HD brings College Park composite expertise to a carbon fiber laminate design. The foot is constructed with aluminum components, and offers a high weight limit. Even with a low-profile and lightweight design, the dynamics of the foot may provide excellent energy return for high-impact K3 users.

HORIZON HD PART NUMBER

HZ SIDE: **L OR R** SHELL COLOR*: **C, T, B, OR J** **EN** SIZE: **21-30** FIRMNESS CATEGORY: **D** **1-7**

* (C)aucasian, (T)an, (B)rown, or (J)et Black. Shell width is standard.
Example part number: HZ LCEN 24 D4

INCLUDED

Endo foot includes: sandal toe foot shell and CPI sock

HORIZON HD FIRMNESS CATEGORIES

PATIENT WEIGHT	0-140lbs	141-180	181-220	221-275	276-330	331-365
SIZE	21-30 cm					26-30 cm
LOW - MODERATE	1	2	3	4	5	7
HIGH IMPACT	2	3	4	5	6	7

TECHNICAL SPECIFICATIONS

FUNCTION	Carbon Fiber Base
SIZES	21-30 cm
MOUNTING	Endo (IP)
WEIGHT LIMIT	330 lbs (21-25 cm) 365 lbs (26-30 cm)
WARRANTY	3 years
BUILD HEIGHT	2.3 in
WEIGHT*	559 g
L-CODE**	IMPACT LEVEL L5981
WATER CAPABILITY Waterproof in Fresh Water	ENVIROHELL*** standard


ENVIROHELL PART NUMBER

HZ S SIDE: **L OR R** SHELL COLOR**: **C, T, B, OR J** SIZE: **21-30**

** (C)aucasian, (T)an, (B)rown, or (J)et Black



*26 cm with Enviroshell. **The listing of codes with these products should not be construed as a guarantee for coverage or payment. Ultimate responsibility for the coding of services/products rests with the individual practitioner. ***Caucasian, Tan, Brown, or Jet Black.



HARDWARE	PART NUMBER	DESCRIPTION
	HZ HWK 2125	Horizon Wedge Kit, 21-25cm
	HZ HWK 2630	Horizon Wedge Kit, 26-30cm
	CP SSET	Super Set Screws (8pc) M8 x 14mm

TOOLS	PART NUMBER	DESCRIPTION
	TFH 800	Foot Horn, Stainless Steel
	CP HK 4	Hex Key, 4mm
	GA GEL	Cyanoacrylate Adhesive Gel (CA glue)

COSMETIC ATTACHMENT PLATE (CAP)	PART NUMBER (LEFT/RIGHT)		DESCRIPTION	ADDITIONAL INFORMATION
	BE CLΔ 2125	BE CRΔ 2125	Cosmetic Attachment Plate (left or right), 21-25cm	Plastic ring snaps into foot shell for easy cosmetic finishing.
	BE CLΔ 2630	BE CRΔ 2630	Cosmetic Attachment Plate (left or right), 26-30cm	

Δ Indicate Color. (C)aucasian, (B)rown or (T)an.

SHELLTREAD	PART NUMBER (LEFT/RIGHT)		DESCRIPTION	ADDITIONAL INFORMATION
	HZ TLN x	HZ TRN x	Shelltread (left or right), Neutral (indicate size)	Increases durability and prevents premature wear of foot shell.

SOCKS	PART NUMBER	DESCRIPTION	ADDITIONAL INFORMATION
	CPI SS SM	CPI Sock, Small, 21-25cm	Highly durable sock that acts as a noise barrier and protects the foot from debris.
	CPI SS LG	CPI Sock, Large, 26-30cm	
	CPI SSK SM	CPI Sock, Small, 21-25cm, 6-pack	
	CPI SSK LG	CPI Sock, Large, 26-30cm, 6-pack	