

ICON[®]

microprocessor knee
 user friendly design
 streamlined setup process
 responses in real time

K3/4



The Icon microprocessor knee offers a user-friendly experience for patients and clinicians. Featuring responsive sensors, streamlined setup, and the intuitive Stride Studio app, this dynamic knee was designed to enhance user confidence and trust. Clinicians can easily access all adjustments, custom modes, and more in their version of Stride Studio. The Icon has no service requirement with College Park's remote updates and unique warranty process. Our engineers have developed patented technology to put the Icon in a league of its own—improving angular accuracy, response time, and weight savings. Equipped with a long-lasting battery and IP68 rating, the Icon is the versatile solution for moderate to high activity users.

PART NUMBER

IC MPKA 1000

TECH SPECS

| | |
|---|--|
| DISTAL ATTACHMENT Pyramid | WEIGHT 1530 g |
| WEIGHT LIMIT 330 lbs 150 kg | FLEXION RANGE 130° |
| OVERALL HEIGHT 10.5 in 267 mm | DOMES TO KNEE CENTER 0.67 in 17 mm |
| WARRANTY 3 years | L-CODES* L5828, L5845, L5848, L5850, L5856, L5925 |

WATER CAPABILITY







FRESH WATER PROOF

*The listing of codes with these products should not be construed as a guarantee for coverage or payment. Ultimate responsibility for the coding of services/products rests with the individual practitioner.


RECOMMENDED ADAPTERS

| | PART NUMBER | DESCRIPTION |
|--|-------------|---|
| | TA T250 30 | 250 mm Low-Profile Tube Adapter, 30 mm diameter, Titanium |
| | TA T430 30 | 430 mm Low-Profile Tube Adapter, 30 mm diameter, Titanium |
| | TC T 30 | Tube Clamp, 30 mm diameter, Titanium |
| | 3PA R | 3-Prong anchor, rotatable |
| | LPTC T | Low Profile Threaded Connector, Ti |

REPLACEMENT PARTS

| PART NUMBER | DESCRIPTION |
|---|---|
| IC WAR3 | Icon Microprocessor Knee 3 Year Extended Warranty |
|  | IC KPK Icon Knee Pad Replacement Kit |
|  | IC BATT Icon Knee Replacement Battery |
|  | IC CHG G Icon Knee Charging Kit, Global |
|  | IC CB G Icon Knee Charging Blade Kit, Global |

HARDWARE & TOOLS

| PART NUMBER | DESCRIPTION |
|--|---------------------------------|
| CA GELK | Knee Cyanoacrylate Adhesive Gel |
|  | CP HK 4 Hex Key, 4mm |