

INTELLIWEAVE COMPOSITE | CURVED HYDRAULIC ANKLE | DUAL-STAGE STRIDE CONTROL



The Odyssey K2 foot was designed to provide fluid control for superior stability and incredible ground contact. The patented curved hydraulic ankle, Intelliweave composite base, and customizable dual Stride Control offer excellent range of motion in a lightweight, low-profile design.

ODYSSEY K2 PART NUMBER

OK2 SIDE: L OR R SHELL COLOR*: C, T OR B **EN** SIZE: 21-30

*(C)aucasian, (T)an or (B)rown. Example part number: OK2 LCEN 24 D3

FIRMNESS CATEGORY

D 1-4

INCLUDED

Endo foot includes: sandal toe foot shell and CPI sock

ODYSSEY K2 FIRMNESS CATEGORIES

PATIENT WEIGHT	0-140lbs	141-180	181-220	221-300
SIZE	21-30 cm			25-30 cm
LOW IMPACT	1	2	3	4

TECHNICAL SPECIFICATIONS

FUNCTION	Curved Hydraulic
SIZES	21-30 cm
MOUNTING	Endo (IP)
WEIGHT LIMIT	220 lbs (21-24 cm) 300 lbs (25-30 cm)
WARRANTY	2 years
BUILD HEIGHT	3.1 in (21-24 cm) 3.3 in (25-30 cm)
WEIGHT*	744 g
L-CODE**	L5972, L5968
WATER CAPABILITY	ENVIROSHELL***
Weatherproof	C T B standard


ENVIROSHELL PART NUMBER





CS S SIDE: L OR R SHELL COLOR**: C, T, OR B SIZE: 21-30


** (C)aucasian, (T)an or (B)rown



*26 cm with Enviroshell. **The listing of codes with these products should not be construed as a guarantee for coverage or payment. Ultimate responsibility for the coding of services/products rests with the individual practitioner. ***Caucasian, Tan, or Brown.



HARDWARE	PART NUMBER	DESCRIPTION
	CP SSET	Super Set Screws (8pc) M8 x 14mm
	TFH 800	Foot Horn, Stainless Steel
	CP HK 3	Hex Key, 3mm, for adjusting Stride Control. Included with foot order.
	CP HK 4	Hex Key, 4mm

SHELLTREAD	PART NUMBER (LEFT/RIGHT)	DESCRIPTION	ADDITIONAL INFORMATION
	CS TLN x CS TRN x	Shelltread (left or right), Neutral (indicate size)	Increases durability and prevents premature wear of foot shell.

SOCKS	PART NUMBER	DESCRIPTION	ADDITIONAL INFORMATION
	CPI SS SM	CPI Sock, Small, 21-25cm	Highly durable sock that acts as a noise barrier and protects the foot from debris.
	CPI SS LG	CPI Sock, Large, 26-30cm	
	CPI SSK SM	CPI Sock, Small, 21-25cm, 6-pack	
	CPI SSK LG	CPI Sock, Large, 26-30cm, 6-pack	