

INTEGRATED SPRING TECHNOLOGY | INTELLIWEAVE COMPOSITE | 10 MM OF VERTICAL SHOCK | WATERPROOF



The first to incorporate Integrated Spring Technology (iST®), the Soleus foot provides natural movement with a progressive, 3-staged response and excellent range of motion. The foot integrates Intelliweave composite springs through a tough, flexible plantar belt that responds to varying impacts. The unique, scaled design may provide 10 mm of vertical compliance for comfort. The Soleus can accommodate a moderate activity lifestyle up to the highest-impact Paralympic athlete.

INCLUDED

Endo foot includes: foot shell, cosmetic attachment plate (CAP), CPI sock and 4-hole fastener kit

SOLEUS PART NUMBER

SO	SIDE	SHELL COLOR*	HOUSING	SIZE	WIDTH**	D	FIRMNESS CATEGORY
	L OR R	C, T, B OR J	M OR T	21-30	N OR W		11-76 (SEE NEXT PAGE)

*(C)aucasian, (T)an, (B)rown or (J)et Black. †Midnight Black or Tactical. **(N)arrow or (W)ide. Wide only available in 24-30cm. Additional cost for Tactical. Example part number: SO LCM 24N D32

ENVIROHELL PART NUMBER

V S	SIDE	SHELL COLOR**	SIZE	WIDTH†
	L OR R	C, T, B OR J	21-30	N OR W


** (C)aucasian, (T)an, (B)rown or (J)et Black. †Wide shell only available in 24-30cm.



TECHNICAL SPECIFICATIONS




FUNCTION	iST®
SIZES	21-30 cm
MOUNTING	Endo (4-hole)
WEIGHT LIMIT	220 lbs (21-24 cm) 250 lbs (25-26 cm) 275 lbs (27-30 cm)
WARRANTY	3 years
BUILD HEIGHT	6.4 in (21-24 cm) 7.1 in (25-26 cm) 7.6 in (27-30 cm)
WEIGHT*	667 g
L-CODE** L5981 PDAC APPROVED	IMPACT LEVEL
WATER CAPABILITY Waterproof	ENVIROHELL*** narrow & wide

*26 cm with Enviroshell. **The listing of codes with these products should not be construed as a guarantee for coverage or payment. Ultimate responsibility for the coding of services/products rests with the individual practitioner. ***Caucasian, Tan, Brown, or Jet Black.




SOLEUS FIRMNESS CATEGORIES

PATIENT WEIGHT	0-100lbs	101-120	121-140	141-160	161-180	181-200	201-220	221-250	251-275
21-24 cm									
MODERATE IMPACT	11	12	33	34	54	55	55		
HIGH IMPACT	12	33	34	54	55	55			
25-26 cm									
MODERATE IMPACT	12	12	32	33	54	55	75	76	
HIGH IMPACT	12	32	33	54	55	75	76	76	
27-30 cm									
MODERATE IMPACT	11	11	12	32	33	34	54	55	55
HIGH IMPACT	12	12	32	33	34	54	55	55	55


HARDWARE	PART NUMBER	DESCRIPTION
	CP SSET	Super Set Screws (8pc) M8 x 14mm
	TFH 800	Foot Horn, Stainless Steel
	CP HK 4	Hex Key, 4mm

RECOMMENDED TITANIUM 4-HOLE ADAPTERS	PART NUMBER	DESCRIPTION
	4HPR T	4-hole Pyramid Receiver, Titanium
	4HPRR T	4-hole Pyramid Receiver, Rotatable, Titanium
	4HPAT	4-hole Pyramid Adapter, Titanium
	4HPAR T	4-hole Pyramid Adapter, Rotatable, Titanium
	4HPH T	4-hole Pyramid Adapter with pin-lock hole, Titanium

RECOMMENDED TITANIUM 4-HOLE ADAPTERS	PART NUMBER	DESCRIPTION
	4HPHR T	4-hole Pyramid Adapter with pin-lock hole, Rotatable, Titanium
	4HPA DOT	4-hole Pyramid Adapter, Diagonal Offset, Titanium
	4HPA AOT	4-hole Pyramid Adapter, Axial Offset, Titanium
	4HFK 14	4-hole Fastener Kit (4pc) M6 x 14mm, Non-Rotatable Adapters
	4HFK 16	4-hole Fastener Kit (4pc) M6 x 16mm, Rotatable Adapters

COSMETIC ATTACHMENT PLATE (CAP)	PART NUMBER (LEFT/RIGHT)	DESCRIPTION	ADDITIONAL INFORMATION
	V CLΔ 2124 V CRΔ 2124	Cosmetic Attachment Plate (left or right), 21-24cm	Plastic ring snaps into foot shell for easy cosmetic finishing.
	V CLΔ 2526 V CRΔ 2526	Cosmetic Attachment Plate (left or right), 25-26cm	
	V CLΔ 2730 V CRΔ 2730	Cosmetic Attachment Plate (left or right), 27-30cm	

Δ Indicate Color. (C)aucasian, (B)rown or (T)an.

SHELLTREAD	PART NUMBER (LEFT/RIGHT)	DESCRIPTION	ADDITIONAL INFORMATION
	V TLN xx V TRN xx	Shelltread (left or right) Neutral (indicate size, width)	Increases durability and prevents premature wear of foot shell.
	V TLJ xx V TRJ xx	Shelltread (left or right) Jet Black (indicate size, width)	

SOCKS	PART NUMBER	DESCRIPTION	ADDITIONAL INFORMATION
	SO SS 2126	CPI Sock, 21-26cm	Highly durable sock that acts as a noise barrier and protects the foot from debris.
	SO SS 2730	CPI Sock, 27-30cm	
	SO SSK 2126	CPI Sock, 6-pack, 21-26cm	
	SO SSK 2730	CPI Sock, 6-pack, 27-30cm	

INSPIRED BY HEROES BUILT FOR HEROES THE TACTICAL SOLEUS

The Tactical Soleus was inspired by Dale Beatty, Military Veteran, bilateral amputee and co-founder of Purple Heart Homes. Driven by the belief that no Veteran should be left behind, his organization provides personalized housing solutions for service-disabled Veterans and their families.



College Park donates \$200 of each Tactical Soleus purchase directly back to Purple Heart Homes. For more information visit phhusa.org

