

INTELLIWEAVE COMPOSITE | SERVICE-FREE DESIGN | CUSTOM GAIT MATCH



The Tempo provides high-end College Park gait matching in a value class foot. The service-free design comes with a full length toe lever, helping to provide stored energy and ground contact for stability. The Tempo is custom-built for each individual and is an ideal solution for K2 patients that need extra durability and controlled range of motion.

TEMPO PART NUMBER

TO **L** OR **R** **C, T** OR **B** **EN** **21-30** **D** **1-5**

*(C)aucasian, (T)an or (B)rown. Shell width is standard. Example part number: TO LCEN 24 D3

INCLUDED




Endo foot includes: bonded sandal toe foot shell


TEMPO FIRMNESS CATEGORIES



PATIENT WEIGHT	0-140lbs	141-180	181-220	221-300	301-330
SIZE	21-30cm				25-30cm
LOW IMPACT	1	2	3	4	5


*25 cm with Enviroshell. **The listing of codes with these products should not be construed as a guarantee for coverage or payment. Ultimate responsibility for the coding of services/products rests with the individual practitioner. ***Caucasian, Tan, and Brown.

TECHNICAL SPECIFICATIONS

FUNCTION	K2
SIZES	21-30 cm
MOUNTING	Endo (IP)
WEIGHT LIMIT	300 lbs (21-24 cm) 330 lbs (25-30 cm)
WARRANTY	2 years
BUILD HEIGHT	2.2 in (21-24 cm) 2.5 in (25-26 cm) 2.6 in (27-30 cm)
WEIGHT*	500 g
L-CODE**	IMPACT LEVEL L5972 
WATER CAPABILITY	ENVIROSHELL***  Weatherproof  standard



HARDWARE	PART NUMBER	DESCRIPTION
	CP SSET	Super Set Screws (8pc) M8 x 14mm
	CP HK 4	Hex Key, 4mm

SHELLTREAD	PART NUMBER (LEFT/RIGHT)	DESCRIPTION	ADDITIONAL INFORMATION
	CS TLN x CS TRN x	Shelltread (left or right), Neutral (indicate size)	Increases durability and prevents premature wear of foot shell.