



metro liners®



COLLEGE PARK
TECHNOLOGY
— FOR THE —
HUMAN RACE
SINCE 1988

FIT AND FLEXIBILITY



The College Park Genesee liner provides the ideal benefits required of a gel liner in an affordable package. Shaped for fit and flexibility with cushioning comfort and excellent limb conforming properties, this durable liner accommodates a wide range of users. The unique fabric provides a superior fiber to gel bond, reduces pistoning and eases donning. The liner is infused with mineral oil to help absorb shear forces and aid in donning and doffing.




MINERAL OIL INFUSED

Protects against shear and friction, as well as aids in donning and doffing.



RETRO-FITS TO INDUSTRY

Retro-fits with industry standard liners of like thickness and profile.



LONGER LENGTH

The Genesee is 16 in (40.6 cm) long for the perfect fit. (+/- .5 in, 1.2 cm)

The Genesee liner features a unique fabric to gel bond, making it more durable than other gel liners.



TECHNICAL SPECS

Material	Thermoplastic Elastomer (TPE) Gel
Types	Cushion or Locking
Sizes	M, M+, L, L+, XL
Activity Level	K1-K3
Length	16 in (+/- .5 cm) 40.6 cm (+/- 1.2 cm)

GEL SLEEVE


The College Park's Genesee Gel Sleeve is suitable as the primary suspension method, or for use with a suction or vacuum suspension. The sleeve is made with hypo-allergenic material and infused with mineral oil to help protect against shear and friction. With a 3mm dimension, the sleeve provides comfort and knee flexion.



With a durable fabric-to-gel bond, the sleeve is able to be trimmed to length without the need for adhesive.




Reinforced interior for puncture-resistance.



Available in 5 sizes and neutral tones.

DISTAL ROLL-ONS

The Genesee Gel Distal Roll-On is suitable to be used alone, or with other gel and fabric-faced prosthetic interfaces. It's made with hypo-allergenic material and infused with mineral oil. The roll-on will self-adhere to other Genesee Gel products.



3 available sizes, 2mm thick, and 8.5" in length to trim to size.



Able to be cut into a band of any width for added comfort under a suspension sleeve.



COLLEGE PARK METRO LINERS SANITIZING SPRAY

Kills 99.9% of all harmful bacteria to diminish odors and sanitize the liner. Available in 4oz, 17oz and 5-pack 17oz bottles.



College Park designs and manufactures a full line of anatomically correct, custom built prosthetic solutions. College Park is ISO 13485 certified with an extensive focus on quality testing, and proudly sources and produces all products entirely in the USA. With precise engineering and quality manufacturing, College Park is committed to innovating human locomotion solutions for users all over the world.



college park
technology for the human race

800.728.7950 | 586.294.7950 | www.college-park.com



facebook.com/CollegeParkIndustries



youtube.com/CollegeParkInd



linkedin.com/company/College-Park-Industries