



College Park designs and manufactures a full line of anatomically correct, custom built prosthetic solutions. College Park is ISO 13485 certified with an extensive focus on quality testing, and proudly sources and produces all products entirely in the USA. With precise engineering and quality manufacturing, College Park is committed to innovating human locomotion solutions for users all over the world.



college park

TECHNOLOGY *for the* HUMAN RACE

800.728.7950 | 586.294.7950 | www.college-park.com



483 LIT BRO SK 200211

©2020, College Park Industries, Inc. All rights reserved.



 sidekicks®



Introducing the Sidekicks™, the first adjustable stubbie feet for foreshortened prostheses. College Park engineers applied 15 years of knowledge gained from the Venture® ankle technology to develop the Sidekicks. With natural anatomical ankle motion and a flexible platform for ground compliance, they are a vast improvement over current options for bilateral, transfemoral amputees. They also come gait-matched with a toolbox of solutions for in-office customization.



WATER FRIENDLY

Sidekicks were tested for water use and are safe to use in the shower or at the beach.



UNEVEN GROUND COMPLIANCE

Flat panel stubbie feet don't accommodate uneven terrain. The Sidekicks were designed to flex and plant for ground contact in real-world environments, crucial for balance and stability.



MOLDED TREAD

Adhered to the foot with a molecular bond, the molded tread on the Sidekicks provide traction and added durability.



CUSTOM GAIT MATCHING

Sidekicks were designed to promote balance and mobility. Swapping bumpers and utilizing the Stride Control™ feature aid in increased ability throughout the rehabilitation process. This is especially important when the goal is a transition to microprocessor knees.



ANKLE ROTATION

The feet mimic anatomical rotation which allow for movement to occur naturally, reducing socket forces and eliminating unwanted torque.



PATIENT TESTIMONIAL

“While competing in New Zealand, I ran them through four airports, three cities, and two ski resorts, through sand, snow, water, mud, and more. They held up to all conditions.”



TECHNICAL SPECIFICATIONS

FUNCTION

Stubbie

MOUNTING

Endo (IP)

CLEARANCE

2.3 in
(5.8 cm)

WEIGHT LIMIT

275 lbs (125 kg)

WEIGHT

422g

WARRANTY

2 years

L CODES*

L5972
L5986

IMPACT LEVEL



*The listing of codes with these products should not be construed as a guarantee for coverage or payment. Ultimate responsibility for the coding of services/products rests with the individual practitioner.