



College Park designs and manufactures a full line of anatomically correct, custom built prosthetic solutions. College Park is ISO 13485 certified with an extensive focus on quality testing, and proudly sources and produces all products entirely in the USA. With precise engineering and quality manufacturing, College Park is committed to innovating human locomotion solutions for users all over the world.



college park

TECHNOLOGY *for the* HUMAN RACE

800.728.7950 | 586.294.7950 | www.college-park.com



496 LIT BRO TP 210419

©2021, College Park Industries, Inc. All rights reserved.



 truper®



The rugged Truper foot combines high energy with stability for active kids. The foot's stable foreheel and energy-storing bumpers provide a smooth, controlled transition from standing to running. Made to easily adjust for growing children, the Truper offers two size ranges that provide flexibility through adolescence.

UNPARALLELED INCLINE/DECLINE



INCLINE: FRONT BUMPER

The front bumper compresses to dorsiflex the ankle and drop the heel for natural movement.



DECLINE: REAR BUMPER

The rear bumper compresses on decline to draw the toe down for additional ground contact and stability.



DESIGNED WITH PLAY IN MIND

Highest grade aluminum provides fatigue resistance.

Integrated bumpers are made of aircraft-grade polymer for longevity.



Specially designed toes with cold flow properties in the material provide curl resistance to prevent an "elfing" effect.

Ridges in the foreheel act like I-beams in a bridge, providing additional strength and reinforcement.

The foreheel was redesigned for increased durability. High-strength nylon fights deformation while providing stability.

GROWTH CAPABILITIES



The Truper can adapt quickly to growing kids with custom growth plate kits and interchangeable foreheels.



TECHNICAL SPECIFICATIONS

FUNCTION

Tri-Axial

MOUNTING

Endo (IP)
Exo (ALX)

SIZES

16-21 cm

WEIGHT LIMIT

16-18 cm	100 lbs (45kg)
19-21 cm	132 lbs (60 kg)

CLEARANCE

2.1-2.4 in (IP) (5.3-6.1 cm)
4.3-4.6 in (ALX) 11.0-11.8 cm

WEIGHT*

215-314g

WARRANTY

1 year

L CODES**

L5979

IMPACT LEVEL



ENVIROSHELL®†



standard

*16-21 cm model, including Enviroshell. **The listing of codes with these products should not be construed as a guarantee for coverage or payment. Ultimate responsibility for the coding of services/products rests with the individual practitioner.

† Caucasian, Tan or Brown.

