

VERTICAL SHOCK ABSORPTION | HIGH-PERFORMANCE COMFORT | MATERIAL FUSION


Soar through life with the Aviator. This lightweight foot is designed to provide high energy return, comfort, and vertical shock absorption in every step. Combining flexible Intelliweave[®] composites with dynamic carbon fiber, the sleek Aviator uses a hybrid material approach in its springs. These features offer superior range of motion to accommodate a wide range of activities and users. When life gets turbulent, the Aviator can help you stay the course.

AVIATOR PART NUMBER

AV	SIDE	SHELL COLOR*	EN	SIZE	D	FIRMNESS CATEGORY
	L OR R	C, T, B OR J		21-30		1-5

*(C)aucasian, (T)an, (B)rown or (J)et Black. ‡(N)arrow or (W)ide. Example part number: AV LBEN 21 D1




INCLUDED

Endo foot includes: foot shell, CPI sock and wedge kit

AVIATOR FIRMNESS CATEGORIES

PATIENT WEIGHT	0-140lbs	141-180	181-220	221-275	276-330
SIZE	21-25cm				
MODERATE-HIGH IMPACT	1	2	3	4	
SIZE	26-30cm				
MODERATE-HIGH IMPACT	1	2	3	4	5


TECHNICAL SPECIFICATIONS

FUNCTION	Vertical Shock
SIZES	21-30 cm
MOUNTING	Endo (IP)
WEIGHT LIMIT	275 lbs (21-25 cm) 330 lbs (26-30 cm)
WARRANTY	3 years
BUILD HEIGHT	6.5 in
WEIGHT*	979 g
L-CODE** L5987 PDAC APPROVED	IMPACT LEVEL 
WATER CAPABILITY  Waterproof in Fresh Water	ENVIROHELL***  standard


ENVIROHELL PART NUMBER



HZ	S	SIDE	SHELL COLOR**	SIZE
		L OR R	C, T, B, OR J	21-30




** (C)aucasian, (T)an, (B)rown, or (J)et Black




*26 cm with Enviroshell. **The listing of codes with these products should not be construed as a guarantee for coverage or payment. Ultimate responsibility for the coding of services/products rests with the individual practitioner. ***Caucasian, Tan, Brown and Jet Black.



HARDWARE	PART NUMBER	DESCRIPTION
	VY HWK 2125	Voyager Wedge Kit, 21-25
	VY HWK 2630	Voyager Wedge Kit, 26-30cm
	CP SSET	Super Set Screws (8pc) M8 x 14mm

TOOLS	PART NUMBER	DESCRIPTION
	TFH 800	Foot Horn, Stainless Steel
	CP HK 4	Hex Key, 4mm
	GA GEL	Cyanoacrylate Adhesive Gel (CA glue)

SHELLTREAD	PART NUMBER (LEFT/RIGHT)	DESCRIPTION	ADDITIONAL INFORMATION
	HZ TLN xx HZ TRN xx	Shelltread (left or right) Neutral (indicate size, width)	Increases durability and prevents premature wear of foot shell.
	HZ TLJ xx HZ TRJ xx	Shelltread (left or right) Jet Black (indicate size, width)	

SOCKS	PART NUMBER	DESCRIPTION	ADDITIONAL INFORMATION
	CPI SS SM	CPI Sock, Small, 21-25cm	Highly durable sock that acts as a noise barrier and protects the foot from debris.
	CPI SS LG	CPI Sock, Large, 26-30cm	
	CPI SSK SM	CPI Sock, Small, 21-25cm, 6-pack	
	CPI SSK LG	CPI Sock, Large, 26-30cm, 6-pack	