

ORBIT®

adjustable torsion
integrated energy return
carbon fiber base



The Orbit foot was designed to provide smooth rotation and rollover. College Park's unique, built-in torsion feature has integrated energy return, which may improve comfort for golf, baseball, tennis, and more! The durable Orbit offers an incredibly low profile and total rotation range of 40°. This foot also incorporates College Park's composite expertise in a carbon fiber design, helping to provide a dynamic response during low to high impact activities.



TECH SPECS

FUNCTION

Torsion

BUILD HEIGHT

3.7 in (21-25 cm)

3.9 in (26-30 cm)

SIZES

21-30 cm

WEIGHT LIMIT

330 lbs

L-CODE†

L5981 + L5984

MOUNTING

Endo (IP)

WEIGHT*

895 g

IMPACT LVL



WARRANTY

3 years

ENVIROSHELL†



standard

*26 cm model, including Enviroshell. †The listing of codes with these products should not be construed as a guarantee for coverage or payment. Ultimate responsibility for the coding of services/products rests with the individual practitioner.

† Caucasian, Tan, or Brown.

ORBIT PART NUMBER

	SIDE	SHELL COLOR*		SIZE		FIRMNESS CATEGORY
OR	L OR R	C, T, B, OR J	EN	21-30	D	1-6

*(C)aucasian, (T)an, (B)rown, or (J)et Black. Shell width is standard.

Example part number: OR LTEN 21 D1

ORBIT FIRMNESS CATEGORIES

PATIENT WEIGHT	0-140lbs	141-180	181-220	221-275	276-330
SIZE	21-30 cm				
LOW-MODERATE IMPACT	1	2	3	4	5
HIGH IMPACT	2	3	4	5	6

INCLUDED

Endo foot includes:
sandal toe foot shell,
CPI sock, and
Cosmetic Attachment
Plate (CAP)


ENVIROSHELL PART NUMBER

	SIDE	SHELL COLOR**	SIZE
HZ S	L OR R	C, T, B, OR J	21-30

** (C)aucasian, (T)an, (B)rown, or (J)et Black.





Enviroshell included with every foot order



HARDWARE AND TOOLS

	PART NUMBER	DESCRIPTION
	HZ HWK 2125	Horizon Wedge Kit, 21-25cm
	HZ HWK 2630	Horizon Wedge Kit, 26-30cm


	PART NUMBER	DESCRIPTION
	CP SSET	Super Set Screws (8pc) M8 x 14mm
	TFH 800	Foot Horn, Stainless Steel

COSMETIC ATTACHMENT PLATE (CAP)

	PART NUMBER (LEFT/RIGHT)	DESCRIPTION	ADDITIONAL INFORMATION
	BE CLΔ 2125 BE CRΔ 2125	Cosmetic Attachment Plate (left or right), 21-25cm	Plastic ring snaps into foot shell for easy cosmetic finishing.
	BE CLΔ 2630 BE CRΔ 2630	Cosmetic Attachment Plate (left or right), 26-30cm	

Δ Indicate Color. (C)aucasian, (B)rown or (T)an.

SHELLTREAD

	PART NUMBER (LEFT/RIGHT)	DESCRIPTION	ADDITIONAL INFORMATION
	HZ TLN x HZ TRN x	Shelltread (left or right), Neutral (indicate size)	Increases durability and prevents premature wear of foot shell.

SOCKS

	PART NUMBER	DESCRIPTION	ADDITIONAL INFORMATION
	HZ SS 2125	CPI Sock, 21-25cm	Highly durable sock that acts as a noise barrier and protects the foot from debris.
	HZ SS 2630	CPI Sock, 26-30cm	
	HZ SSK 2125	CPI Sock, 6-pack, 21-25cm	
	HZ SSK 2630	CPI Sock, 6-pack, 26-30cm	