



Skinergy Plus

**19**  
**TrueFinish™**  
**Colour Shades**



## Skinergy Leg Covers

Lifelike, durable and easy to fit, SkinergyPlus offers a fantastic, off-the-shelf, cosmetic finish for the below knee amputee.

Available in 19 skin shades and featuring a NEW silicone formulation for a more natural finish.



# Skinergy Plus

skinergyplus

## Features:

- Integrated reinforced mesh
- 19 colour shades
- Seamless design
- Easy to clean
- Available in closed or split toe
- Adult foot sizes 22-30 cm
- Two calf sizes

| Split Toe<br>Part Number | Closed Toe<br>Part Number | Foot Size<br>(cm) | Ankle Circ. |      | Mid Calf Circ. |      | Top of Cover<br>Circ. |      | Cover Length From Floor to - |           |
|--------------------------|---------------------------|-------------------|-------------|------|----------------|------|-----------------------|------|------------------------------|-----------|
|                          |                           |                   | Min         | Max  | Min            | Max  | Min                   | Max  | Mid Calf                     | Cover Top |
| TLC22LS/E*               | TNS22LS/E*                | 22-23             | 20          | 25   | 28             | 35   | 26                    | 32.5 | 32                           | 55        |
| TLC22RS/E*               | TNS22RS/E*                | 22-23             | 20          | 25   | 28             | 35   | 26                    | 32.5 | 32                           | 55        |
| TLC22LL/E*               | TNS22LL/E*                | 22-23             | 20          | 25   | 32             | 40   | 30                    | 37.5 | 32                           | 55        |
| TLC22RL/E*               | TNS22RL/E*                | 22-23             | 20          | 25   | 32             | 40   | 30                    | 37.5 | 32                           | 55        |
| TLC24LS/E*               | TNS24LS/E*                | 24-25             | 22          | 27.5 | 31             | 39   | 29                    | 36   | 34.5                         | 57        |
| TLC24RS/E*               | TNS24RS/E*                | 24-25             | 22          | 27.5 | 31             | 39   | 29                    | 36   | 34.5                         | 57        |
| TLC24LL/E*               | TNS24LL/E*                | 24-25             | 22          | 27.5 | 35             | 44   | 33                    | 41   | 34.5                         | 57        |
| TLC24RL/E*               | TNS24RL/E*                | 24-25             | 22          | 27.5 | 35             | 44   | 33                    | 41   | 34.5                         | 57        |
| TLC26LS/E*               | TNS26LS/E*                | 26-27             | 22          | 28   | 32             | 40   | 30                    | 37.5 | 37                           | 61        |
| TLC26RS/E*               | TNS26RS/E*                | 26-27             | 22          | 28   | 32             | 40   | 30                    | 37.5 | 37                           | 61        |
| TLC26LL/E*               | TNS26LL/E*                | 26-27             | 22          | 28   | 37             | 46   | 35                    | 44   | 37                           | 61        |
| TLC26RL/E*               | TNS26RL/E*                | 26-27             | 22          | 28   | 37             | 46   | 35                    | 44   | 37                           | 61        |
| TLC28LS/E*               | TNS28LS/E*                | 28-30             | 24          | 30   | 34             | 42.5 | 32                    | 40   | 39.5                         | 66        |
| TLC28RS/E*               | TNS28RS/E*                | 28-30             | 24          | 30   | 34             | 42.5 | 32                    | 40   | 39.5                         | 66        |
| TLC28LL/E*               | TNS28LL/E*                | 28-30             | 24          | 30   | 39             | 49   | 37                    | 46   | 39.5                         | 66        |
| TLC28RL/E*               | TNS28RL/E*                | 28-30             | 24          | 30   | 29             | 49   | 37                    | 46   | 39.5                         | 66        |

## 19 TrueFinish™ Colour Shades



| Part Number | Description                             |
|-------------|-----------------------------------------|
| NCT9740     | Skinergy Plus TrueFinish 25mm Punchtool |
| E24655      | Skinergy Plus TrueFinish Swatch Kits    |
| E26625/E*   | Skinergy Plus TrueFinish Patch Kits     |
| E26626/E*   | Skinergy Plus TrueFinish Spare Parts    |

**Note:** E\* denotes colour. For example, TNS22LS/E5