












# ENDOSKELETAL COMPONENTS REFERENCE

T 800.728.7950 586.294.7950 | F 800.294.0067 586.294.0067 | college-park.com



## TUBE ADAPTERS



COMPONENT	BUILD HEIGHT	WEIGHT	WEIGHT LIMIT	PART NUMBER
 <b>WATERPROOF TUBE ADAPTER, Titanium</b> 250 mm 30 mm diameter	95 mm (min)	255 g (max)	300 lbs	TA WT250 30
 <b>WATERPROOF TUBE ADAPTER, Titanium</b> 430 mm 30 mm diameter	95 mm (min)	390 g (max)	300 lbs	TA WT430 30
 <b>LOW-PROFILE TUBE ADAPTER, Titanium</b> 250 mm 30 mm diameter	57 mm (min)	221 g (max)	300 lbs	TA T250 30
 <b>LOW-PROFILE TUBE ADAPTER, Titanium</b> 430 mm 30 mm diameter	57 mm (min)	337 g (max)	300 lbs	TA T430 30
 <b>LOW-PROFILE TUBE ADAPTER, Stainless Steel</b> 250 mm 30 mm diameter	57 mm (min)	303 g (max)	300 lbs	TA S250 30
 <b>LOW-PROFILE TUBE ADAPTER, Stainless Steel</b> 430 mm 30 mm diameter	57 mm (min)	418 g (max)	300 lbs	TA S430 30

COMPONENT	BUILD HEIGHT	WEIGHT	WEIGHT LIMIT	PART NUMBER
 <b>LOW-PROFILE TUBE ADAPTER, Stainless Steel</b> 250 mm 34 mm diameter	57 mm (min)	302 g (max)	400 lbs	TA S250 34
 <b>LOW-PROFILE TUBE ADAPTER, Stainless Steel</b> 430 mm 34 mm diameter	57 mm (min)	425 g (max)	400 lbs	TA S430 34
 <b>LOW-PROFILE TUBE ADAPTER, Titanium</b> 300 mm 25 mm diameter with 30 mm adult receiver	52 mm (min)	242 g (max)	165 lbs	TA T300 2530
 <b>LOW-PROFILE TUBE ADAPTER, Titanium</b> 300 mm 25 mm diameter	47 mm (min)	209 g (max)	165 lbs	TA T300 25
 <b>LOW-PROFILE TUBE ADAPTER, Titanium</b> 300 mm 22 mm diameter	47 mm (min)	169 g (max)	132 lbs	TA T300 22




## TUBE CLAMPS




COMPONENT	BUILD HEIGHT	WEIGHT	WEIGHT LIMIT	PART NUMBER
 <b>TUBE CLAMP, Titanium</b> 30 mm diameter	51 mm	82 g	300 lbs	TC T 30
 <b>TUBE CLAMP, Stainless Steel</b> 30 mm diameter	51 mm	146 g	300 lbs	TC S 30
 <b>TUBE CLAMP, Stainless Steel</b> 34 mm diameter	55 mm	302 g	400 lbs	TC S 34
 <b>TUBE CLAMP, Titanium</b> 25 mm diameter with 30 mm adult receiver	45 mm	83 g	165 lbs	TC T 2530
 <b>TUBE CLAMP, Titanium</b> 25 mm diameter	39 mm	45 g	165 lbs	TC T 25
 <b>TUBE CLAMP, Titanium</b> 22 mm diameter	39 mm	46 g	132 lbs	TC T 22

## DOUBLE ENDED ADAPTERS



COMPONENT	BUILD HEIGHT	WEIGHT	WEIGHT LIMIT	PART NUMBER
 <b>DOUBLE ENDED ADAPTER, Titanium</b> 35 mm	35 mm	95 g	300 lbs	DR T35
 <b>DOUBLE ENDED ADAPTER, Titanium</b> 45 mm	45 mm	110 g	300 lbs	DR T45
 <b>DOUBLE ENDED ADAPTER, Titanium</b> 60 mm	60 mm	123 g	300 lbs	DR T60
 <b>DOUBLE ENDED ADAPTER, Titanium</b> 75 mm	75 mm	136 g	300 lbs	DR T75




## TORSION ADAPTERS

COMPONENT	BUILD HEIGHT	WEIGHT	WEIGHT LIMIT	PART NUMBER
 <b>TORSION ADAPTER, PYRAMID ADAPTER TO PYRAMID RECEIVER</b> Titanium	43 mm	393 g	330 lbs	TOR PAPT
 <b>TORSION ADAPTER, TUBE CLAMP TO PYRAMID RECEIVER,</b> Stainless Steel	47 mm	432 g	330 lbs	TOR TCPT
 <b>TORSION ADAPTER, 30MM TUBE ADAPTER, TI, 275MM,</b> Titanium (receiver), Aluminum (tube)	83 mm (min)	567 g	330 lbs	TOR T275 30

COMPONENT	BUILD HEIGHT	WEIGHT	WEIGHT LIMIT	PART NUMBER
 <b>TORSION ADAPTER, 30MM TUBE ADAPTER, TI, 455MM,</b> Titanium (receiver), Aluminum (tube)	83 mm (min)	670 g	330 lbs	TOR T455 30
 <b>TORSION ADAPTER, 34MM TUBE ADAPTER, TI, 275MM,</b> Titanium (receiver), Aluminum (tube)	79 mm (min)	569 g	330 lbs	TOR T275 34
 <b>TORSION ADAPTER, 34MM TUBE ADAPTER, TI, 455MM,</b> Titanium (receiver), Aluminum (tube)	79 mm (min)	685 g	330 lbs	TOR T455 34

## FLEXION PLATES

COMPONENT	BUILD HEIGHT	WEIGHT	WEIGHT LIMIT	PART NUMBER
 <b>4-HOLE THREADED, 10° W/ ADAPTER W/ PINLOCK HOLE, ROTATABLE</b> Stainless Steel & Aluminum	13.4 mm	215 g	300 lbs	4HTP 10 AHR S
 <b>4-HOLE THREADED, 10° W/ RECEIVER, ROTATABLE,</b> Stainless Steel & Aluminum	23.4 mm	214 g	300 lbs	4HTP 10 RR S








COMPONENT	BUILD HEIGHT	WEIGHT	WEIGHT LIMIT	PART NUMBER
 <b>4-HOLE THREADED, 10° W/ ADAPTER W/ PINLOCK HOLE, ROTATABLE,</b> Titanium & Aluminum	13.4 mm	187 g	300 lbs	4HTP 10 AHR T
 <b>4-HOLE THREADED, 10° W/ RECEIVER, ROTATABLE,</b> Titanium & Aluminum	23.4 mm	186 g	300 lbs	4HTP 10 RR T
 <b>4-HOLE THREADED, 10°,</b> Aluminum	10.2 mm	149 g	300 lbs	4HTP 10 A







# ENDOSKELETAL COMPONENTS REFERENCE

T 800.728.7950 586.294.7950 | F 800.294.0067 586.294.0067 | college-park.com













## 4-HOLE PYRAMID ADAPTERS AND RECEIVERS

COMPONENT	BUILD HEIGHT	WEIGHT	WEIGHT LIMIT	PART NUMBER
 4-HOLE PYRAMID RECEIVER, Titanium Standard	17 mm	64 g	300 lbs	4HPR T
 4-HOLE PYRAMID RECEIVER, Stainless Steel Standard	17 mm	75 g	400 lbs	4HPR S
 4-HOLE PYRAMID RECEIVER, Titanium Standard	15 mm	44 g	165 lbs	4HPR TP
 4-HOLE PYRAMID RECEIVER, Titanium Rotatable	19 mm	90 g	300 lbs	4HPRR T
 4-HOLE PYRAMID RECEIVER, Titanium Rotatable	16 mm	50 g	165 lbs	4HPRR TP
 4-HOLE PYRAMID ADAPTER, Titanium Standard	7 mm	46 g	300 lbs	4HPA T
 4-HOLE PYRAMID ADAPTER, Stainless Steel Standard	7 mm	118 g	400 lbs	4HPA S





COMPONENT	BUILD HEIGHT	WEIGHT	WEIGHT LIMIT	PART NUMBER
 4-HOLE PYRAMID ADAPTER, Titanium Standard	8 mm	36 g	165 lbs	4HPA TP
 4-HOLE PYRAMID ADAPTER, Titanium Rotatable	9.5 mm	77 g	300 lbs	4HPAR T
 4-HOLE PYRAMID ADAPTER, Titanium Rotatable	9 mm	56 g	165 lbs	4HPAR TP
 4-HOLE PYRAMID OFFSET, Titanium Axial	7 mm	55 g	300 lbs	4HPA AO T
 4-HOLE PYRAMID OFFSET, Titanium Diagonal	7 mm	64 g	300 lbs	4HPA DO T
 4-HOLE PYRAMID ADAPTER, Titanium w/Pin-Lock Hole	7 mm	44 g	300 lbs	4HPH T
 4-HOLE PYRAMID ADAPTER, Titanium w/Pin-Lock Hole, Rotatable	7 mm	74 g	300 lbs	4HPHR T





## LAMINATION ADAPTERS

COMPONENT	BUILD HEIGHT	WEIGHT	WEIGHT LIMIT	PART NUMBER
 4-EAR SOCKET ADAPTER W/ PIN LOCK HOLE, Titanium	5.7 mm	55.2 g	300 lbs	4ESAH T
 4-EAR SOCKET ADAPTER W/ PIN LOCK HOLE, Stainless Steel	5.7 mm	96.5 g	300 lbs	4ESAH S
 THREADED ADAPTER W/PIN LOCK HOLE, Titanium	3.2 mm	38 g	300 lbs	TAH T
 THREADED ADAPTER W/PIN LOCK HOLE, Stainless Steel	3.2 mm	66.4 g	300 lbs	TAH S
 THREADED RECEIVER, Titanium	13.2 mm	36.8 g	300 lbs	TR T
 THREADED RECEIVER, Stainless Steel	13.2 mm	64.5 g	300 lbs	TR S
 LOW PROFILE THREADED CONNECTOR, Titanium	n/a	33 g	330 lbs	LPTC T
 3-PRONG ANCHOR, ROTATABLE, Stainless Steel	8 mm	124 g	365 lbs	3PA R
 4-PRONG ANCHOR, ROTATABLE, Stainless Steel	8 mm	158 g	365 lbs	4PA R

COMPONENT	BUILD HEIGHT	WEIGHT	WEIGHT LIMIT	PART NUMBER
 4-HOLE PLATE, SQUARE, THREADED, Aluminum	6.4 mm	36 g	300 lbs	4HLP ST A
 4-HOLE PLATE, ROUND, THREADED, TIE GROOVE, Aluminum	9.7 mm	59 g	300 lbs	4HLP RTG A
 4-HOLE PLATE, ROUND, THREADED, TIE GROOVE, Aluminum	9.7 mm	62 g	300 lbs	4HLP RG A
 4-HOLE PLATE, ROUND, THREADED, Aluminum	6.4 mm	54 g	300 lbs	4HLP RT A
 4-HOLE PLATE, SQUARE, Aluminum	6.4 mm	40 g	300 lbs	4HLP A
 4-HOLE PLATE, ROUND, HEAVY DUTY, Titanium	6.4 mm	36 g	300 lbs	4HLP RHD T
 4-HOLE PLATE, ROUND, HEAVY DUTY, Stainless Steel	6.4 mm	68 g	300 lbs	4HLP RHD S
 4-HOLE PLATE, SQUARE, PEDIATRIC, Aluminum	6.4 mm	31 g	300 lbs	4HLP AP

## LAMINATION KITS

COMPONENT	BUILD HEIGHT	WEIGHT	WEIGHT LIMIT	PART NUMBER
 3 PRONG ANCHOR, ADAPTER W/ PIN LOCK HOLE, ROTATABLE, Titanium	11.2 mm	162 g	300 lbs	3PA AHR T
 3 PRONG ANCHOR W/ RECEIVER, ROTATABLE, Titanium	21.2 mm	160 g	300 lbs	3PA RR T
 4 PRONG ANCHOR, ADAPTER W/ PIN LOCK HOLE, ROTATABLE, Titanium	11.2 mm	196 g	300 lbs	4PA AHR T
 4 PRONG ANCHOR W/ RECEIVER, ROTATABLE, Titanium	21.2 mm	195 g	300 lbs	4PA RR T

COMPONENT	BUILD HEIGHT	WEIGHT	WEIGHT LIMIT	PART NUMBER
 3 PRONG ANCHOR, ADAPTER W/ PIN LOCK HOLE, ROTATABLE, Stainless Steel	11.2 mm	190 g	300 lbs	3PA AHR S
 3 PRONG ANCHOR W/ RECEIVER, ROTATABLE, Stainless Steel	21.2 mm	188 g	300 lbs	3PA RR S
 4 PRONG ANCHOR, ADAPTER W/ PIN LOCK HOLE, ROTATABLE, Stainless Steel	11.1 mm	224 g	300 lbs	4PA AHR S
 4 PRONG ANCHOR W/ RECEIVER, ROTATABLE, Stainless Steel	21.2 mm	223 g	365 lbs	4PA RR S