

# KNEES



# GUARDIAN™

## friction brake knee

easy adjustments  
stance flexion  
remote lock



The Guardian knee is a lightweight mechanical knee joint featuring a friction brake for stance control, stance flexion, and extension assist. The knee is best suited as a rehabilitation tool for new patients to learn standing and walking during the rehabilitation process, and as a permanent knee for low-impact patients. It provides the prosthetist easy adjustability for patient gait-matching without the need to remove the prosthesis. The Guardian also has a remote lock feature that can be disengaged as the patient progresses to provide the full 145 degrees of anatomical motion.

### PART NUMBER

# GFKA 120

### TECH SPECS

#### MATERIAL

Aluminum

#### WEIGHT

537 g

#### WEIGHT LIMIT

275 lbs

#### FLEXION RANGE

145°

#### WARRANTY

2 years

#### L-CODES\*

L5812, L5845,  
L5850, L5925

\*The listing of codes with these products should not be construed as a guarantee for coverage or payment. Ultimate responsibility for the coding of services/products rests with the individual practitioner.

### HARDWARE & TOOLS

	PART NUMBER	DESCRIPTION
	CP KELK	Knee External Lock Kit (included with knee)
	CP LK	Lanyard Kit (included with knee)
	CP HK 4	Hex Key, 4mm

# VICTOR™

## pneumatic knee

easy adjustments  
stance flexion  
variable cadence



Step victoriously. The Victor knee is a single axis, pneumatic knee joint featuring a friction brake to control stance phase. This knee is the ideal solution for low to moderate activity users, as it can adjust to accommodate both slow and fast paced walking. It also works well for patients looking to increase their activity level after rehabilitation. The Victor is a safe, economical option that instills confidence and freedom in patients.

### PART NUMBER



# VPKA 200

### HARDWARE & TOOLS

### TECH SPECS

<b>MATERIAL</b>	<b>WEIGHT</b>
Aluminum	710 g
<b>WEIGHT LIMIT</b>	<b>FLEXION RANGE</b>
275 lbs	145°
<b>WARRANTY</b>	<b>L-CODES*</b>
2 years	L5822, L5845, L5850

\*The listing of codes with these products should not be construed as a guarantee for coverage or payment. Ultimate responsibility for the coding of services/products rests with the individual practitioner.

	PART NUMBER	DESCRIPTION
	CP KSB	Knee Stop Bumper Replacement Kit
	CP HK 4	Hex Key, 4mm