

		PART NUMBER						FIRMNESS CATEGORY		
FOOT	PART ID	SIDE	SHELL COLOR	MOUNTING		SIZE	WIDE OPTION (24-30cm)	(REFER TO INDIVIDUAL SELECTION CHARTS)		
							WIDE OPTION (23-26cm)			
K3 FEET	ACCENT	AC	L / R	C / B	IP / DF / DM / D4	M / L	21-28	N / W	D	1-3
	AVIATOR	AV	L / R	C / T / B / J	EN		21-30		D	1-5
	HORIZON HD	HZ	L / R	C / T / B / J	EN		21-30		D	1-7
	ODYSSEY iQ	OD iQ	L / R	C / T / B / J	EN		23-30		D	1-5
	ODYSSEY K3	ODK3	L / R	C / T / B / J	EN		21-30		D	1-5
	ORBIT	OR	L / R	C / T / B / J	EN		21-30		D	1-6
	SCOUT (adapter required)	SC	L/R	C / T / B	M		13-22		D	1-4
					HOUSING COLOR					
	SOLEUS	SO	L / R	C / T / B / J	M / T		21-30	N / W	D	11-76
	TERRAIN	TN	L / R	C / T / B / J	EN		21-30		D	1-7
	TRIBUTE (regular toe)	TB	L / R	C / T / B / J	EN		21-30	N / W	D	1-4
	TRIBUTE (sandal toe)	TB	L / R	C / T / B	EN	S	21-30		D	1-4
	TRUPER	TP	L / R	C / T / B	EN / ALX		16-21		D	1-3
	TRUSTEP	TS	L / R	C / T / B / J	EN / ALX		22-31	N / W	D	1-9
	K2 FEET	VELOCITY	VL	L / R	C / T / B / J	M	EN	25-30	N / W	D
VENTURE		V	L / R	C / T / B / J	EN / ALX		21-30	N / W	D	1-5
VOYAGER		VY	L / R	C / T / B / J	EN		21-30		D	1-5
BREEZE		BE	L / R	C / T / B	EN		21-30		D	3
CELSUS		CS	L / R	C / T / B	EN		21-30		D	1-4
CELSUS iVT (tube kit required)		CS iVT	L / R	C / T / B	EN		21-30		D	1-4
ODYSSEY K2		OK2	L / R	C / T / B	EN		21-30		D	1-4
			SINGLE= 1 PAIR= 2							
SIDEKICKS		SK	AK	1/2			ONE SIZE		D	1-3
TEMPO		TO	L / R	C / T / B	EN		21-30		D	1-5
TEMPO iVT (tube kit required)		TO iVT	L / R	C / T / B	EN		21-30		D	1-5
TERRAIN iW		TN iW	L / R	C / T / B	EN		21-30		D	1-4

PRODUCT KEY

HOUSING COLOR		SHELL COLOR		MOUNTING OPTION		BUTTON PLACEMENT	
Midnight Black	M	Caucasian	C	Endo	EN	Medial	M
Tactical	T	Tan	T	Exo	ALX	Lateral	L
		Brown	B	Integrated Pyramid	IP		
		Jet Black	J	DYNAMIC PYLON + ADAPTER			
				DP + Female	DF		
				DP + Male	DM		
				DP + 4-Hole	D4		

K3 FEET - FIRMNESS CATEGORIES

ACCENT

WEIGHT LBS	0-120	121-140	141-160	161-180	181-220	221-250
SIZE CM	21-28					25-28
LOW IMPACT	1	1	2	2	3	3
MODERATE IMPACT	1	2	2	3	3	3

AVIATOR

WEIGHT LBS	0-140	141-180	181-220	221-275	276-330
SIZE CM	21-25				
MODERATE-HIGH IMPACT	1	2	3	4	N/A
SIZE CM	26-30				
MODERATE-HIGH IMPACT	1	2	3	4	5

HORIZON HD

WEIGHT LBS	0-140	141-180	181-220	221-275	276-330	331-365
SIZE CM	21-30					26-30
LOW-MODERATE IMPACT	1	2	3	4	5	7
HIGH IMPACT	2	3	4	5	6	7

ODYSSEY iQ

WEIGHT LBS	0-140	141-180	181-220	221-275
SIZE CM	23-30			
LOW IMPACT	1	2	3	4
MODERATE IMPACT	2	3	4	5

ODYSSEY K3

WEIGHT LBS	0-140	141-180	181-220	221-275
SIZE CM	21-30			
LOW IMPACT	1	2	3	4
MODERATE IMPACT	2	3	4	5

ORBIT

WEIGHT LBS	0-140	141-180	181-220	221-275	276-330
SIZE CM	21-30				
LOW-MODERATE IMPACT	1	2	3	4	5
HIGH IMPACT	2	3	4	5	6

K3 FEET - FIRMNESS CATEGORIES (CONTINUED)

SCOUT

WEIGHT LBS	0-55	56-74	75-100	101-132
SIZE CM	13-15			
HIGH IMPACT	1	N/A	N/A	N/A
SIZE CM	16-18			
HIGH IMPACT	1	2	3	N/A
SIZE CM	19-22			
HIGH IMPACT	1	2	3	4

REQUIRED ADAPTER

(SCOUT ONLY)

Threaded Adapter, SS, Pediatric	TA SP
Threaded Receiver, SS, Pediatric	TR SP
Threaded Exo Adapter, SS, Pediatric**	TXA SP
Threaded Lamination Adapter, Small, SS, Pediatric	TLAS SP
Threaded Lamination Adapter, Medium, SS, Pediatric	TLAM SP
Threaded Lamination Adapter, Large, SS, Pediatric	TLAL SP

Order adapter as a separate line item.

SCOUT FOOT SHELL SIZE GUIDE

FOOT SIZE (CM)	13	14	15	16	17	18	19	20	21	22
Regular Shell (cm)	13	14	15	16	17	18	19	20	21	22
Size UP Shell (cm)	13/14	14/15	N/A	16/17	17/18	N/A	19/20	20/21	21/22	N/A

Check out College Park's unique size-up foot shells!

Contact Customer Service for assistance with ordering Size Up foot shells.

SOLEUS

WEIGHT LBS	0-100	101-120	121-140	141-160	161-180	181-200	201-220	221-250	251-275
SIZE CM	21-24								
MODERATE IMPACT	11	12	33	34	54	55	55	N/A	N/A
HIGH IMPACT	12	33	34	54	55	55	N/A	N/A	N/A
SIZE CM	25-26								
MODERATE IMPACT	12	12	32	33	54	55	75	76	N/A
HIGH IMPACT	12	32	33	54	55	75	76	76	N/A
SIZE CM	27-30								
MODERATE IMPACT	11	11	12	32	33	34	54	55	55
HIGH IMPACT	12	12	32	33	34	54	55	55	55

TERRAIN

WEIGHT LBS	0-140	141-180	181-220	221-275	276-330	331-365
SIZE CM	21-30					26-30
LOW-MODERATE IMPACT	1	2	3	4	5	7
HIGH IMPACT	2	3	4	5	6	7

TRIBUTE

WEIGHT LBS	0-100	101-120	121-140	141-160	161-180	181-200	201-220
SIZE CM	21-30						
LOW IMPACT	1	1	2	2	3	3	4
MODERATE IMPACT	1	2	2	2	3	4	4

K3 FEET - FIRMNESS CATEGORIES (CONTINUED)

TRUPER

WEIGHT LBS	0-49	50-74	75-100	101-132
SIZE CM	16-21			19-21
HIGH IMPACT	1	2	3	3

TRUSTEP

WEIGHT LBS	0-74	75-100	101-120	121-140	141-160	161-180	181-200	201-220	221-250	251-275	276-300	301-330	331-352
SIZE CM	22-31											EXO ONLY	
LOW IMPACT	1	2	3	4	4	5	6	6	7	7	8	8	9
MODERATE IMPACT	1	2	3	4	5	6	6	7	7	7	8	9	9
HIGH IMPACT	2	3	4	5	6	6	7	7	8	8	8	9	9

VELOCITY

WEIGHT LBS	0-100	101-120	121-140	141-160	161-180	181-200	201-220	221-250	251-275
SIZE CM	25-26								
LOW-MODERATE IMPACT	1	1	1	2	2	3	4	5	N/A
HIGH IMPACT	1	2	2	3	3	4	5	6	N/A
SIZE CM	27-30								
LOW-MODERATE IMPACT	1	1	1	1	2	2	3	4	5
HIGH IMPACT	1	2	2	2	3	3	4	5	6

VENTURE

WEIGHT LBS	0-100	101-120	121-140	141-160	161-180	181-200	201-220	221-250	251-275
SIZE CM	21-24								
LOW IMPACT	1	1	2	3	3	4	4	5	N/A
MODERATE IMPACT	1	2	3	3	4	4	5	5	N/A
HIGH IMPACT	2	3	3	3	4	4	5	5	N/A
SIZE CM	25-30								
LOW IMPACT	1	1	2	3	3	4	4	5	5
MODERATE IMPACT	1	2	3	3	4	4	5	5	N/A
HIGH IMPACT	2	3	3	4	4	4	5	5	N/A

VOYAGER

WEIGHT LBS	0-140	141-180	181-220	221-275	276-330
SIZE CM	21-25				
MODERATE-HIGH IMPACT	1	2	3	4	N/A
SIZE CM	26-30				
MODERATE-HIGH IMPACT	1	2	3	4	5

K2 FEET - FIRMNESS CATEGORIES

BREEZE

WEIGHT LBS	0-220	221-275
SIZE CM	21-30	26-30
LOW IMPACT	3	3

CELSUS / CELSUS iVT

WEIGHT LBS	0-140	141-180	181-220	221-300
SIZE CM	21-30			25-30
LOW IMPACT	1	2	3	4

REQUIRED TUBE KIT (iVT ONLY)

iVT AK Tube Kit	IVT AK TBK
iVT BK Tube Kit	IVT BK TBK

selection of one tube kit is required as a separate line item.

ODYSSEY K2

WEIGHT LBS	0-140	141-180	181-220	221-300
SIZE CM	21-30			25-30
LOW IMPACT	1	2	3	4

SIDEKICKS

WEIGHT LBS	0-120	121-220	221-275
SIZE CM	ONE SIZE		
LOW IMPACT	1	2	3

TEMPO / TEMPO iVT

WEIGHT LBS	0-140	141-180	181-220	221-300	301-330
SIZE CM	21-30				25-30
LOW IMPACT	1	2	3	4	5

REQUIRED TUBE KIT (iVT ONLY)

iVT AK Tube Kit	IVT AK TBK
iVT BK Tube Kit	IVT BK TBK

selection of one tube kit is required as a separate line item.

TERRAIN iW

WEIGHT LBS	0-140	141-180	181-220	221-300
SIZE CM	21-30			25-30
LOW IMPACT	1	2	3	4

FOOT SHELL PART NUMBERS (Foot shells are included with every foot purchase)

FOOT	FOOT SHELL	SHELL TYPE	PART ID	SHELL	SIDE	SHELL COLOR	SIZE	WIDE OPTION (24-30 cm)
Breeze	BREEZE	Sandal Toe	BE	S	L / R	C / T / B	21-30	n/a
Celsus, Celsus iVT, Odyssey K2, Tribute, Terrain iW	CELSUS	Sandal Toe	CS	S	L / R	C / T / B	21-30	n/a
Aviator, Horizon HD, Odyssey K3, Orbit, Terrain, Voyager	HORIZON	Sandal Toe	HZ	S	L / R	C / T / B / J	21-30	n/a
Odyssey iQ	ODYSSEY iQ	Sandal Toe	OD iQ	S	L / R	C / T / B / J	23-30	n/a
Scout	SCOUT	Sandal Toe	SC	XS* / S	L / R	C / T / B	13-22**	n/a
Truper	TRUPER	Sandal Toe	TP	XS* / S	L / R	C / T / B	16-21	n/a
Trustep	TRUSTEP	Regular Toe	TS	XS* / S	L / R	C / T / B / J	22-31	N/W
Soleus, Venture, Tribute, Velocity	VENTURE	Regular Toe	V	S	L / R	C / T / B / J	21-30†	N/W

*XS = exo; add 1 day for processing.

**Available Size-Up Footshell Sizes are: 1314, 1415, 1617, 1718, 1920, 2021, 2122

†Velocity size is 25-30.