

INSTRUCTION SHEET

ENDOSKELETAL COMPONENT

30 mm Low Profile Tube Adapter, SS, 250 mm / Part No. TA S250 30

TECHNICAL INFORMATION

Material	Stainless Steel, Aluminum
Build Height (MAX)	250 mm (9.8 in)
Cut Length (MIN)	57 mm (2.24 in)
Assembly Weight (MAX)	303 g
Patient Weight Limit	136 kg (300 lbs)
Warranty	3 Year, limited
Torque (Set Screws)	15 N·m (11 ft-lbs)

M8 x 14mm set screws with nylon patch are included with the adapter.
Loctite® not required.

- ⚠ CAUTIONS**
- Torque evenly to specified settings.
 - Do not expose this product to corrosive materials.

For further questions, please contact College Park Technical Service.



college park

800.728.7950 | 586.294.7950 | college-park.com

MADE IN THE USA

454 INS TA S250 30 170804



COLLEGE PARK INDUSTRIES
27955 College Park Dr.
Warren, MI 48088 USA



EMERGO EUROPE
Prinsessegracht 20, 2514 AP
The Hague, Netherlands

Australian Sponsor
EMERGO AUSTRALIA
Level 20, Tower II, Darling Park,
201 Sussex Street,
Sydney, NSW 2000 Australia

©2017 College Park Industries, Inc. All rights reserved. College Park is a registered trademark of College Park Industries, Inc.