

# Battery Adapters & Mounts

## Battery Cable Adapters:

Below is a partial list of the common Battery Adapters for use with LTI Built-in Battery Recharge Modules. These are the most common Battery Adapters for connecting LTI batteries to other manufacturer's systems. For a full list of battery adapters, see the Battery Adapter section.



**MC31** Battery Adapter, Bock Quick Disconnect 2-socket plug

*Use with any of the BPXX Recharge Modules shown above to connect with any Bock-style power plug*



**MC37** Battery Adapter, 1 Bock Electrode, VASI Plug

*Use with any of the BPXX Recharge Modules shown above to connect with any VASI prosthetic hand*



**MC38** Battery Adapter, 2 Bock Electrodes, VASI Plug

*Use with any of the BPXX Recharge Modules shown above to connect with any VASI prosthetic hand*



**MC41** Battery Adapter, Bock QD 2-socket plug, 5V Regulator

*Use with any of the BPXX Recharge Modules shown above to connect with the Otto Bock 2000 series children's prosthetic hand*

## ***Battery Adapters & Mounts***

### ***Removable Battery Mounts:***

Below are the battery receptacles required for the removable batteries shown above. Each of these receptacles are supplied with a cable and Bock-style 2-socket power connector.



**MC72** Battery Receptacle for BP72 Removable Battery  
*Supplied with cable and 2-socket Bock-style power connector*



**MC90** Battery Receptacle for BP90 Removable Battery  
*Supplied with cable and 2-socket Bock-style power connector*