



college park

technology for the human race

800.728.7950 | 586.294.0067



facebook.com/CollegeParkIndustries



instagram.com/CollegeParkInd



youtube.com/CollegeParkInd



college-park.com/accent



Australian Sponsor
EMERGO AUSTRALIA
Level 20, Tower II, Darling Park,
201 Sussex Street
Sydney, NSW 2000, Australia

EC REP **EMERGO EUROPE**
Prinsessegracht 20, 2514 AP, The Hague, Netherlands

COLLEGE PARK INDUSTRIES, INC.
27955 College Park Dr., Warren, MI 48088 USA



337 LIT RK AC J 170921



I AM

a world traveler

peer mentor

home chef

motivational speaker

nature lover

music & art enthusiast

getting my life back.

I AM COLLEGE PARK.



college park

technology for the human race

I AM JEN



...THIS IS MY STORY.

During her Junior year of college in 2001, Jen lost control while riding a minibike, crushing her ankle, leg, foot and big toe. In the decade after the accident, she had 21 surgeries and stubborn, excruciating staph infections. Jen was emotionally and physically exhausted. Her leg was worse than useless - it was destroying her life. After years of debilitating, chronic pain and severe depression, Jen was determined to find another way. She discovered inspiring stories online of active, positive amputees and learned of the technology available. After much thought, she went in for the amputation surgery almost exactly 10 years to the day after the accident. Today, Jen lives a full and inspiring life. She wears a College Park Soleus foot for activities like bike riding, golfing, and hiking and an

Accent foot anytime she wants to be fashionable and wear heels.



“Since my accident, I’ve only been able to wear sneakers or flats. I’d have half a row to choose from in a shoe store and it was one more thing that reminded me of my disability. Boots are so popular and in style, and I was sad I had to settle. I never want to be held back by my disability, least of all not by SHOES! The Accent gives me that freedom back and having a shoe choice makes me feel feminine and confident. I love dressing up and I own it when I walk into a room now!”

IS THE ACCENT RIGHT FOR YOU?



The Accent foot is best for someone with a moderate activity lifestyle that wants to look and feel great every day. The heel can adjust up to 2 inches to accommodate an entire range of footwear, from dress shoes to flat sandals to cowboy boots. It's the perfect foot for those who care most about beautiful cosmetic finishing and shoe style.



FOOT FEATURES

DESIGNED FOR FASHION. ENGINEERED FOR LIFE.

The Accent was built with beauty in mind. An easy-to-push button is located on the ankle to adjust the heel height and lock it into place. The unique design allows for 13 precise settings to accommodate even the smallest shift in heel height.



COSMETIC OPTIONS

The Accent includes the sandal toe feature and soft ankle fairing that locks into the shell to prevent buckling and provide a contoured ankle.



INTELLIWEAVE COMPOSITE

College Park's proprietary Intelliweave composite technology brings the maximum amount of strength and flexibility.



TECH SPECS

GEAR	WEIGHT	RULER	IMPACT	HEIGHT	WEIGHT
MOUNTING	WEIGHT LIMIT	SIZE RANGE	IMPACT LEVEL	CLEARANCE	WEIGHT*
Endo	220 lbs (21-24 cm)	21-28 cm	L M	3.3-3.9 in (IP)	608 g (IP)
Dynamic Pylon	250 lbs (25-28 cm)			7.8-8.5 in (DP)	761 g (DP)

*INCLUDING ENVIRONMENTAL