



college park

technology for the human race

800.728.7950 | 586.294.7950



facebook.com/CollegeParkIndustries



instagram.com/CollegeParkInd



youtube.com/CollegeParkInd



college-park.com/horizon



Australian Sponsor
EMERGO AUSTRALIA
Level 20, Tower II
Darling Park
201 Sussex Street
Sydney, NSW 2000
Australia

COLLEGE PARK INDUSTRIES, INC
27955 College Park Dr., Warren, MI 48088 USA

EMERGO EUROPE
Prinsessegracht 20, 2514 AP, The Hague, Netherlands



482 LIT RK HZ 170731



I AM

a basketball player
this close to dunking
high school graduate
loyal boyfriend
horse enthusiast
mom's hero
figuring it out.

I AM COLLEGE PARK.

 horizon **HD**



college park

technology for the human race


I AM MICHAEL



...THIS IS MY STORY.

Michael was born with fibular hemimelia, and his right leg below the knee was amputated at 6 months old. Fitted in his first prosthesis at 9 months, Michael was walking by 11 months, jumping on a trampoline at 2 years and - with his mother strongly advocating he is "NOT handicapped" - continues to push limits and inspire people every day. He has a passion for playing sports and playing them well, from baseball, basketball, football, roller blading, even skateboarding and Airsoft. His family has

lived on a horse farm in rural Michigan since Michael was 4 years old and he helps out by hauling hay, carrying feed, scooping manure and other farm chores; preferably ones he needs their four-wheeler to accomplish. He even learned to drive a stick-shift at the age of 9 thanks to his dad. Now 19, Michael plans to travel to Japan and Europe and is thinking about pursuing a career in computers, or possibly as an O&P technician as a way to give back and help other amputees. As Michael often says, "They don't give you a leg so you can just sit there".



The Horizon gives me the great energy return that I need every day, but it is more comfortable and has a much smoother roll-over than my other foot. I also love the sandal toe, it's the first foot I've ever had with one"

IS THE HORIZON RIGHT FOR YOU?

horizon[®] HD stay flexible

The HorizonHD brings high energy function in a compact design, perfect for those with long residual limbs. College Park used their composite expertise to maximum the flexibility in the carbon fiber laminate. This provides a high dynamic response during sports, plus a smooth gait and comfort during everyday walking. Ultra-low and ultra-light weight, the Horizon brings a new element to your game.




FOOT FEATURES



DESIGNED TO MOVE

The carbon fiber heel spring and toe spring were meticulously crafted utilizing maximum stress predictions to increase strength and resist failure.



TITANIUM COMPONENTRY

The HorizonHD is built with titanium components, a lighter, stronger alternative to cheaper stainless steel often used in other feet.



SMOOTH GAIT ENVIROSHIELD

The HorizonHD shell's toes were specially designed to assist with rollover for the smoothest possible gait.



ULTRA-LOW PROFILE

With a clearance of 2.3 in (5.8cm), the HorizonHD is one of the lowest profile, high-impact K3 carbon feet on the planet.



TECH SPECS

MOUNTING	WEIGHT LIMIT	SIZE RANGE	IMPACT LEVEL	CLEARANCE	WEIGHT*
Endo	330 lbs (21-24 cm) 365 lbs (25-30 cm)	21-30 cm	L M H	2.3"	518 g

*INCLUDING ENVIROSHIELD