



# college park

technology for the human race

800.728.7950 | 586.294.0067



facebook.com/CollegeParkIndustries



instagram.com/CollegeParkInd



youtube.com/CollegeParkInd



college-park.com/sidekicks



Australian Sponsor  
**EMERGO AUSTRALIA**  
Level 20, Tower II, Darling Park,  
201 Sussex Street  
Sydney, NSW 2000, Australia

EC REP **EMERGO EUROPE**  
Prinsessegracht 20, 2514 AP, The Hague, Netherlands

**COLLEGE PARK INDUSTRIES, INC.**  
27955 College Park Dr., Warren, MI 48088 USA



MADE IN THE U.S.A.  
535 LIT RK SK 170412



a retired Marine Sergeant  
movie buff  
tenacious  
grateful to be alive  
monoskiier  
training for the paralympics  
determined to succeed.

**I AM COLLEGE PARK.**



# college park

technology for the human race



# I AM JOSH



## ...THIS IS MY STORY.

Marine Sergeant Josh Elliott was gravely injured from stepping on a 20-pound explosive device while on routine patrol in Afghanistan in April 2011. Before his accident, Josh was a master swimmer, ran five days per week, taught martial arts, and had a passion for snowboarding. He did not lose consciousness following the blast and remembers everything, including being transported by helicopter, and realizing that his body had been riddled with shrapnel and that both of his legs above the knees and

three fingers on his left hand were gone. As he recovered, Josh began to try some of the things he had once loved and had done so well. Frustrated with his snowboarding attempts, Josh tried monoskiing and immediately loved it. Within a few months, he tried out and was accepted on the Paralympic development team.

Along with his everyday feet, Josh wears his Sidekicks stubbie feet for activities when he wants to be closer to the ground, like yardwork or a quick transfer. As a rehabilitated and very active amputee, Josh tests the limits of his Sidekicks.



*While competing in New Zealand, I ran them through four airports, three cities, and two ski resorts, through sand, snow, water, mud, and more. They held up to all conditions."*



## ARE SIDEKICKS RIGHT FOR YOU?



All Heroes need a Sidekick.

Sidekicks® are the first multi-axial stubbie feet for bilateral, above-knee amputees. With natural anatomical ankle motion and a flexible platform, they are a vast improvement over current options. Gait-matched and approved for water use, Sidekicks flex and plant in real-world environments for the ultimate balance and stability.



### FOOT FEATURES



#### MOLDED TREAD

Adhered to the foot with a molecular bond, the molded tread on the Sidekicks provide traction and added durability.



#### CUSTOM GAIT MATCHING

Your Prosthetist can adjust the Sidekicks throughout the rehabilitation process for increased ability.



#### UNEVEN GROUND COMPLIANCE

The combination of a flexible platform and anatomical ankle motion give the Sidekicks excellent ground compliance.



#### WATER FRIENDLY

Sidekicks were tested for water use and are safe to use in the shower or at the beach.



### TECH SPECS



MOUNTING

Endo



WEIGHT LIMIT

275 lbs



SIZE RANGE

ONE SIZE



IMPACT LEVEL

C L M H



CLEARANCE

2.3 in



WEIGHT

422 g